Package 'espadon'

October 26, 2022

```
Type Package
Title Easy Study of Patient DICOM Data in Oncology
Version 1.1.1
Description Exploitation, processing and 2D-3D visualization of DICOM-
     RT files (structures, dosimetry, imagery) for medical physics and clinical research, in a patient-
     oriented perspective.
License GPL-3
URL https://espadon.cnrs.fr
Encoding UTF-8
Roxygen list(markdown = TRUE)
RoxygenNote 7.2.0
Imports colorspace,
     DT,
     graphics,
     grDevices,
     igraph,
     js,
     mathjaxr,
     Matrix,
     methods,
     misc3d,
     openxlsx,
     progress,
     Rcpp (>= 1.0.7),
     Rdpack,
     rgl (>= 0.107.14),
     Rvcg,
     shiny,
     shinyWidgets,
     sp,
     stats
LinkingTo Rcpp
RdMacros Rdpack,
     mathjaxr
Depends R (>= 4.0.0)
```

Suggests knitr, rmarkdown

VignetteBuilder knitr

R topics documented:

get.line									
get.obj.connection									
get.plane									
get.rigid.M									
$get.roi.connection\ .\ .\ .\ .$									
get.value.from.ijk	 		. 67						
get.value.from.xyz	 		. 68						
get.volume.from.bin	 		. 69						
get.volume.from.roi	 		. 70						
get.xyz.from.index	 		. 71						
grid.equal	 		. 71						
histo.2D	 		. 72						
histo.DVH	 		. 74						
histo.from.bin	 		. 75						
histo.from.roi	 		. 77						
histo.vol									
load.obj.data	 		. 80						
load.obj.from.dicom	 		. 81						
load.obj.from.Rdcm									
load.patient.from.dicom .									. 83
load.patient.from.Rdcm .									. 84
load.Rdcm.raw.data									. 85
load.T.MAT									. 86
mesh.from.bin									. 87
mesh.in.new.ref									. 89
mesh.repair									. 90
mesh.spheric.proj									. 90
nesting.bin									. 92
nesting.cube									. 93
nesting.roi									. 94
orientation.create									. 95
pal.RVV									. 96
Rdcm.inventory									
ref.add									
ref.cutplane.add									
ref.remove									
ref.srctodest.add									
rt.chi.index									
rt.gamma.index									
rt.indices.from.bin									
rt.indices.from.roi									
save.T.MAT									
save.to.Rdcm									
select.names									
struct.clustering									
struct.from.bin									
struct.in.new.ref									
struct.merge									
toy.dicom.raw									
toy.load.patient									
vector.product									
vol.copy	 	٠	. 124						

4 add.margin

vol.create																				125	5
vol.from.bin																				126	5
vol.gradient																				127	7
vol.in.new.ref .																				128	3
vol.median																				128	3
vol.regrid																				129)
vol.sum																				130)
xlsx.from.dcm .																				131	1
xlsx.from.Rdcm																				132	2

Index 134

add.margin	Adding or removing a margin to a volume	

Description

The add.margin function adds or subtracts a margin of the rectangular parallelepiped circumscribed by a volume.

Usage

```
add.margin(vol, xyz.margin, alias = "", description = NULL)
```

Arguments

vol	"volume" class object.
xyz.margin	Vector of the 3 positive or negative x, y and z margins in mm, in the frame of reference of volume cut planes.
alias	Character string, \$alias of the created object
description	Character string, describing the created object. If description = NULL (default value), it will be set to vol\$description

Value

Returns a "volume" class object (see espadon.class for class definitions), in which 3D volume is restricted or increased by the requested margins. If the created volume exceeds the initial volume, new voxels are set to NA.

See Also

nesting.cube

bin.closing 5

bin.closing

Binary volume closing

Description

The bin.closing function performs a morphological operation of closing, using a sphere, on a "volume" class object of "binary" modality. Closing is useful for :

- filling holes that are smaller than the radius,
- merging two shapes close to each other.

Usage

```
bin.closing(vol, radius = 10, alias = "", description = NULL)
```

Arguments

vol	"volume" class object, of "binary" modality
radius	Positive number, in millimeters. By default, radius = 10.
alias	Character string, \$object.alias of the created object.
description	Character string, describing the created object. If description = NULL (default value), it will be set to paste (vol\$object.alias, "closing r =", radius).

Value

Returns a "volume" class object of "binary" modality (see espadon.class for class definitions), with the same grid as vol, in which \$vol3D.data has been transformed by the closing operation.

Note

Closing can be time consuming, try to reduce the binary volume to the strict minimum, before any operations.

See Also

bin.dilation, bin.erosion, bin.opening, add.margin, nesting.cube.

6 bin.clustering

Examples

bin.clustering

Binary volume clustering

Description

The bin.clustering function groups and labels TRUE voxels that have a 6-connectivity (i.e. sharing a common side).

Usage

```
bin.clustering(vol, alias = "", description = NULL)
```

Arguments

vol "volume" class object, of "binary" modality
alias Character string, \$alias of the created object.
description Character string, describing the created object. If description = NULL (default

value), it will be set to paste (vol\$object.alias, "clustering")

Value

Returns a "volume" class object (see espadon.class for class definitions), of "cluster" modality. This object contains the \$cluster.info field, detailing the label and volumes in cm^3 of the different clusters. Note that the label "0" is used for the background.

bin.dilation 7

bin.dilation

Binary volume dilation

Description

The bin.dilation function enlarges a "volume" class object, of "binary" modality, by means of convolution with a sphere. Dilation is useful for :

- filling holes that are smaller than the radius,
- enlarging capes,
- filling narrow channels,
- merging two shapes close to each other.

Usage

```
bin.dilation(vol, radius = 10, alias = "", description = NULL)
```

Arguments

vol	"volume" class object, of "binary" modality
radius	Positive number, in millimeters. By default, radius = 10.
alias	Character string, \$object.alias of the created object.
description	Character string, describing the created object. If description = $NULL$ (default value), it will be set to paste (vol\$object.alias, "dilataion r =", radius)

Value

Returns a "volume" class object of "binary" modality (see espadon.class for class definitions), with the same grid as vol, in which the selected volume has been enlarged by the radius.

Note

Dilation can be time consuming, try to reduce the binary volume to the strict minimum, before any operations.

See Also

bin.erosion, bin.opening, bin.closing, add.margin, nesting.cube.

8 bin.erosion

Examples

bin.erosion

Binary volume erosion

Description

The bin.erosion function decreases a "volume" class object, of "binary" modality, by means of convolution with a sphere. Erosion is useful for :

- removing volumes that are smaller than the radius,
- eliminating narrow capes,
- enlarging channels,
- turning peninsulas into islands.

Usage

```
bin.erosion(vol, radius = 10, alias = "", description = NULL)
```

Arguments

vol	"volume" class object, of "binary" modality
radius	Positive number, in millimeters. By default, radius = 10.
alias	Character string, \$object.alias of the created object.
description	Character string, describing the created object. If description = NULL (default value), it will be set to paste (vol\$object.alias, "erosion r = ", radius).

Value

Returns a "volume" class object of "binary" modality (see espadon.class for class definitions), with the same grid as vol, in which the selected volume has been reduced by the radius.

Note

Erosion can be time consuming, try to reduce the binary volume to the strict minimum, before any operations.

bin.from.roi 9

See Also

bin.dilation, bin.opening, bin.closing, add.margin, nesting.cube.

Examples

bin.from.roi

Creation of a binary volume according to RoI

Description

The bin.from.roi function creates a "volume" class object, of "binary" modality, in which all the voxels of a RoI are set to TRUE.

Usage

```
bin.from.roi(
  vol,
  struct,
  roi.name = NULL,
  roi.sname = NULL,
  roi.idx = NULL,
  T.MAT = NULL,
  within = TRUE,
  alias = "",
  description = NULL
```

Arguments

vol "volume" class object.

struct "struct" class object.

roi.name Vector of exact names of the RoI in the struct object. By default roi.name = NULL. See Details.

roi.sname Vector of names or parts of names of the RoI in the struct object. By default roi.sname = NULL. See Details.

10 bin.from.roi

roi.idx	Vector of indices of the RoI that belong to the struct object. By default roi.idx = NULL. See Details.
T.MAT	"t.mat" class object, created by load.patient.from.Rdcm or load.T.MAT. If T.MAT = NULL, struct\$ref.pseudo must be equal to vol\$ref.pseudo or set to NULL.
within	Boolean, defaults to TRUE. If within = TRUE, the contours included in a RoI are managed, depending on their \$level field. If within = FALSE, only the \$level = 0 fields of the RoI are used (i.e. only the external outlines).
alias	Character string, \$alias of the created object.
description	Character string, describing the created object. If description = NULL (default value), it will be set to struct\$roi.info\$roi.pseudo[roi.idx].

Details

roi.name, roi.sname, and roi.idx must select only one RoI.

Value

Returns a "volume" class object of "binary" modality (see espadon.class for class definitions), with the same grid as vol, in which the voxels in the RoI are set to TRUE.

See Also

bin.from.vol.

```
# loading of toy-patient objects (decrease dxyz for better result)
step <- 3
patient <- toy.load.patient (modality = c("ct", "rtstruct"),</pre>
                              roi.name = c("eye", "optical nerve", "brain"),
                              dxyz = rep (step, 3))
CT <- patient$ct[[1]]
S <- patient$rtstruct[[1]]</pre>
# "optical nerve" binary without inclusions management
bin <- bin.from.roi (CT, struct = S, roi.sname = "left optical",</pre>
                     alias = "left_optical_nerve")
display.plane (CT, top = bin, struct = S,
               view.coord = S$roi.info[S$roi.info$roi.pseudo == "leftopticalnerve",]$Gz,
               legend.shift = -80, interpolate = FALSE, main = "Left nerve selection")
## Not run:
# with a smaller step
patient <- toy.load.patient (modality = c("ct", "rtstruct"),</pre>
                              roi.name = c("eye", "optical nerve", "brain"),
                              dxyz = rep (step, 3))
CT <- patient$ct[[1]]</pre>
S <- patient$rtstruct[[1]]</pre>
# "optical nerve" binary without inclusions management
bin <- bin.from.roi (CT, struct = S, roi.sname = "left optical",</pre>
                      alias = "left_optical_nerve", within = FALSE)
display.plane (CT, top = bin, struct = S,
               view.coord = S$roi.info[S$roi.info$roi.pseudo == "leftopticalnerve",]$Gz,
```

bin.from.vol

bin.from.vol

Creation of a binary volume according to the voxel values of a volume

Description

The bin.from.vol function creates a "volume" class object, of "binary" modality, in which the voxels fulfilling a condition on their value are selected.

Usage

```
bin.from.vol(
  vol,
  min = -Inf,
  max = Inf,
  in.selection = TRUE,
  alias = "",
  description = NULL
)
```

Arguments

vol "volume" class object.
min Minimum value of the selected voxel. Default to -Inf.
max Maximum value of the selected voxel. Default to +Inf.
in.selection Boolean, defaults to TRUE. If in.selection = FALSE, the selected pixels are those whose value is not between min and max.
alias Character string, \$alias of the created object.
description Character string, describing the created object. If description = NULL (default value), it will be set to paste (min, vol\$object.alias, max, sep = "<="") or if in.selection = FALSE, paste ("!(", description, ")").</pre>

Value

Returns a "volume" class object of "binary" modality, with the same grid as vol, in which the selected voxels (i.e. set to TRUE) are those satisfying the following conditions:

- If in.selection = TRUE, then min <= vol\$vol3D.data <= max.
- If in.selection = FALSE, then vol\$vol3D.data < min or max < vol\$vol3D.data

12 bin.intersection

Examples

bin.intersection

Intersection of two binaries

Description

The bin.intersection function creates a "volume" class object, of "binary" modality, representing the intersection (logical AND) of two binary objects.

Usage

```
bin.intersection(vol1, vol2, alias = "", description = NULL)
```

Arguments

```
vol1, vol2 "volume" class objects, of "binary" modality.

alias Character string, $alias of the created object.

description Character string, describing the created object. If description = NULL (default value), it will be set to paste (vol1$object.alias, "&", vol2$object.alias).
```

Value

Returns a "volume" class object of "binary" modality (see espadon.class for class definitions), with the same grid as vol1 and vol2, intersection of vol1 and vol2.

bin.inversion 13

bin.inversion

Inversion of a binary

Description

The bin.inversion function creates a "volume" class object, of "binary" modality, representing the inverse (logical NOT) of another binary object.

Usage

```
bin.inversion(vol, alias = "", description = NULL)
```

Arguments

vol "volume" class object, of "binary" modality

alias Character string, \$alias of the created object.

description Character string, describing the created object. If description = NULL (default value), it will be set to paste ("!", vol\$object.alias, sep = "").

Value

Returns a "volume" class object of "binary" modality (see espadon.class for class definitions), with the same grid as vol, inverse of vol.

14 bin.opening

k	oin.opening	Binary volume opening	

Description

The bin opening function performs a morphological operation of opening, using a sphere, on a "volume" class object of "binary" modality. Opening is useful for:

- removing volumes that are smaller than the radius,
- · smoothing shapes.

Usage

```
bin.opening(vol, radius = 10, alias = "", description = NULL)
```

Arguments

vol	"volume" class object, of "binary" modality.
radius	Positive number, in millimeters. By default, radius = 10.
alias	Character string, \$object.alias of the created object.
description	Character string, describing the created object. If description = $NULL$ (default value), it will be set to paste (vol\$object.alias, "opening r =", radius).

Value

Returns a "volume" class object of "binary" modality (see espadon.class for class definitions), with the same grid as vol, in which \$vol3D.data has been transformed by the opening operation.

Note

Opening can be time consuming, try to reduce the binary volume to the strict minimum, before any operations.

See Also

bin.dilation, bin.erosion, bin.closing, add.margin, nesting.cube.

bin.subtraction 15

bin.subtraction

Subtraction of two binaries

Description

The bin. subtraction function creates a "volume" class object of "binary" modality, representing the subtraction of two binary objects.

Usage

```
bin.subtraction(vol1, vol2, alias = "", description = NULL)
```

Arguments

```
vol1, vol2 "volume" class objects of "binary" modality.

alias Character string, $alias of the created object.

description Character string, describing the created object. If description = NULL (default value), it will be set to paste (vol1$object.alias, "-", vol2$object.alias).
```

Value

Returns a "volume" class object of "binary" modality (see espadon.class for class definitions), with the same grid as vol1 and vol2, in which vol2 is subtracted from vol1.

```
# loading of toy-patient objects (decrease dxyz for better result)
step <- 4
patient <- toy.load.patient (modality = c("mr", "rtstruct"), roi.name = "",</pre>
                               dxyz = rep (step, 3))
MR <- patient$mr[[1]]</pre>
S <- patient$rtstruct[[1]]</pre>
z.ptv <- S$roi.info$Gz[S$roi.info$roi.pseudo == "ptv"]</pre>
# binaries
bin.patient <- bin.from.roi (MR, struct = S, roi.name = "patient",</pre>
                            alias = "patient", T.MAT = patient$T.MAT)
bin.ptv <- bin.from.roi (MR, struct = S, roi.name = "ptv",</pre>
                            alias = "ptv", T.MAT = patient$T.MAT)
#' calculation of the 'patient - ptv' binary
bin <- bin.subtraction (bin.patient, bin.ptv, alias = "patient - ptv")</pre>
display.plane (MR, top = bin, view.coord = z.ptv,
                display.ref = S$ref.pseudo, T.MAT = patient$T.MAT,
                interpolate = FALSE)
```

16 bin.sum

bin.sum

Sum of two binaries

Description

The bin.sum function creates a "volume" class object of "binary" modality, representing the sum (logical OR) of two binary objects.

Usage

```
bin.sum(vol1, vol2, alias = "", description = NULL)
```

Arguments

Value

Returns a "volume" class object of "binary" modality (see espadon.class for class definitions), with the same grid as vol1 and vol2, sum of vol1 and vol2.

```
# loading of toy-patient objects (decrease dxyz for better result)
step <- 4
patient <- toy.load.patient (modality = c("ct", "rtstruct"), roi.name = "eye",</pre>
                               dxyz = rep (step, 3))
CT <- patient$ct[[1]]</pre>
S <- patient$rtstruct[[1]]</pre>
z.leye<- S$roi.info$Gz[S$roi.info$roi.pseudo == "lefteye"]</pre>
# 'left eye' et 'right eye' binaries
bin.left.eye <- bin.from.roi (CT, struct = S, roi.sname = "lefteye",</pre>
                                alias = "left eye")
bin.right.eye <- bin.from.roi (CT, struct = S, roi.name = "righteye",</pre>
                                 alias = "right eye")
bin.eyes <- bin.sum (bin.left.eye, bin.right.eye, alias = "eyes")</pre>
display.plane (CT, top = bin.eyes, struct = S, roi.sname = "eye",
                view.coord = z.leye, legend.shift = -90 ,
                interpolate = FALSE)
```

castlow.str 17

castlow.str

Cast of a character string

Description

The castlow.str function converts a word to lowercase, without accents and spaces.

Usage

```
castlow.str(st)
```

Arguments

st

character string

Value

Returns the ASCII//TRANSLIT transcription of the word st, without accents, spaces and in lower-case letters.

See Also

castup.str.

Examples

```
castlow.str (st = c("Right eye", "Left_Lung", "Right-Lung"))
```

castup.str

Cast of a character string

Description

The castup.str function converts a word to upper case, without accents and spaces.

Usage

```
castup.str(st)
```

Arguments

st

character string

Value

Returns the ASCII//TRANSLIT transcription of the word st, without accents, spaces and in capitals.

See Also

castlow.str.

18 dicom.browser

Examples

```
castup.str (st = c("Right eye", "Left_Lung", "Right-Lung"))
```

dicom.browser

DICOM raw data browser

Description

the dicom.browser function creates a dataframe describing the tags contained in the raw data of a DICOM file, as well as the information to access them.

Usage

```
dicom.browser(
  dicom.raw.data,
  nbTAG = 0,
  stop.tag = "",
  stop.level = 0,
  full.info = FALSE,
  tag.dictionary = dicom.tag.dictionary()
)
```

Arguments

```
dicom.raw.data Raw vector, representing the binary extraction of the DICOM file.

nbTAG Integer. If nbTAG = 0 (default), and stop.tag = "", all the DICOM raw data is browsed. Otherwise, the function only browses the first nbTAG tags.

stop.tag Character string, representing the tag that stops the browse of the dicom.raw.data.

stop.level Positive integer, specifying the encapsulation level of the stop.tag in dicom.raw.data.

full.info Boolean. If TRUE, more information about the DICOM data is returned.

tag.dictionary Dataframe, by default equal to dicom.tag.dictionary, whose structure it must keep. This dataframe is used to parse DICOM files.
```

Value

Returns a dataframe if dicom.raw.data is DICOM raw data, NULL otherwise.

If full.info = FALSE, dataframe columns are

- tag: the tags contained in dicom.raw.data,
- VR : value representation of the content of the tag,
- endian: the endianness of the tag content,
- start: the start address in dicom.raw.data of the tag content.
- stop: the stop address in dicom. raw. data of the tag content.

If full.info = TRUE, the following columns are added:

- encaps.load: If the tag contains nested data, this column gives the number of bytes remaining until the end of the nesting. If there are several levels of nesting, these numbers are collapsed and separated by a space.
- load.start: the start address in dicom.raw.data of the tag load size.
- load.stop: the stop address in dicom.raw.dataof the tag load size.
- tag.start: the start address in dicom.raw.data of the tag.

dicom.parser 19

See Also

dicom.raw.data.loader, dicom.tag.parser

Examples

```
\# DICOM information dataframe of the dummy raw data toy.dicom.raw () df <- dicom.browser (toy.dicom.raw (), full.info = TRUE) str (df)
```

dicom.parser

Conversion of DICOM raw data into a dataframe or a list of DICOM TAG information

Description

The dicom.parser function creates a dataframe or a list from DICOM raw data. The created dataframe or list provides information about the content of the DICOM TAGs included in the raw data.

Usage

```
dicom.parser(
  dicom.raw.data,
  as.txt = TRUE,
  nested.list = FALSE,
  try.parse = FALSE,
  txt.sep = "\\",
  txt.length = 100,
  tag.dictionary = dicom.tag.dictionary()
)
```

Arguments

dicom.raw.data	Raw vector, representing the binary extraction of the DICOM file.
as.txt	Boolean. If as.txt = TRUE, the function returns a dataframe, a list otherwise.
nested.list	Boolean. Only used if as.txt = FALSE. If nested.list = FALSE, the returned list consists of nested lists.
try.parse	Boolean. If TRUE, the tag with unknown DICOM \mbox{VR} (value representation) is converted into string if possible.
txt.sep	String. Used if as.txt = TRUE. Separator of the tag value elements.
txt.length	Positive integer. Used if as.txt = TRUE. Maximum number of letters in the representation of the TAG value.
tag.dictionary	Dataframe, by default equal to dicom.tag.dictionary, whose structure it must keep. This dataframe is used to parse DICOM files.

Value

Returns a list of elements or a dataframe, depending on as.list.

If it returns a dataframe, the columns are names TAG, VR (value representation), VM (value multiplicity), loadsize and Value. The field \$Value is a string representation of the true value.

If it returns a list, each of its elements, named by a TAG, is either a vector or a string, depending of the TAG included in dicom.raw.data.

See Also

dicom.raw.data.loader, dicom.tag.parser

Examples

```
# content of the dummy raw data toy.dicom.raw (), as a list.
L <- dicom.parser (toy.dicom.raw (), as.txt = FALSE)
str(L[40:57])

L <- dicom.parser (toy.dicom.raw (), as.txt = FALSE, nested.list = TRUE)
str(L[40:45])

# content of the dummy raw data toy.dicom.raw (), as a dataframe.
L <- dicom.parser (toy.dicom.raw (), as.txt = TRUE)
str (L)</pre>
```

dicom.raw.data.anonymizer

DICOM anonymizer

Description

the dicom.raw.data.anonymizer function anonymizes dicom.raw.data.

Usage

```
dicom.raw.data.anonymizer(
   dicom.raw.data,
   offset = 0,
   new.PIN = "Anonymous ",
   reset.private.tag = FALSE,
   tag.dictionary = dicom.tag.dictionary()
)
```

Arguments

dicom.raw.data Raw vector, representing the binary extraction of the DICOM file.

offset Integer, default to 0. Each date of the DICOM will be shifted by this offset

expressed in days.

new.PIN Character string, representing the PIN remplacing the old one.

reset.private.tag

Boolean, if TRUE, the value of tags that are not in the tag. $\operatorname{dictionary}$ is removed.

dicom.raw.data.loader 21

tag.dictionary Dataframe, by default equal to dicom.tag.dictionary, whose structure it must keep. This dataframe is used to parse DICOM files.

Value

Returns an anonymyzed raw vector. See Note.

Note

The raw data is anonymized as follows:

- Each date of the DICOM file will be shifted by offset expressed in days.
- Each patient's name, and patient'ID are remplaced by new.PIN
- · All other patient data are deleted, except age, weight, height, gender and shifted birthday.
- All address, phone, physician, operator, author, reviewer, service.
- If reset.private.tag = TRUE, the values of the tags not contained in the tag.dictionary are deleted.

Examples

dicom.raw.data.loader DICOM file loading in raw data

Description

the dicom.raw.data.loader function loads a DICOM file as raw data.

Usage

```
dicom.raw.data.loader(dcm.filename)
```

Arguments

dcm. filename Character string, representing the full name of a DICOM file.

22 dicom.set.tag.value

Value

Returns a vector of raw data from dcm. filename.

See Also

dicom.browser, dicom.tag.parser

Examples

dicom.set.tag.value Change TAG value in DICOM raw data

Description

The dicom.set.tag.value function changes, in the DICOM raw data, the values of the TAG whose VR is a string of characters.

Usage

```
dicom.set.tag.value(
  dicom.raw.data,
  tag,
  tag.value,
  tag.dictionary = dicom.tag.dictionary()
)
```

Arguments

dicom.raw.data Raw vector, representing the binary extraction of the DICOM file.

tag String vector, representing the list of tags whose value is to be changed. See

note 1.

tag.value String vector,representing the list of new tag values.

tag.dictionary Dataframe, by default equal to dicom.tag.dictionary, whose structure it must

keep. This dataframe is used to parse DICOM files.

dicom.tag.dictionary 23

Value

Returns a raw vector, with new tag values.

Note

1- The list of tags included in the DICOM file are given by the first columns of the dataframe provided by the functions dicom.browser and dicom.parser.

2- The dicom.set.tag.value function may take some processing time. To minimize this time, it is recommended to prepare in advance all the tags to be modified, and use the dicom.set.tag.value function only once, as shown in the example.

Examples

dicom.tag.dictionary DICOM TAG dictionary

Description

The dicom.tag.dictionary function gives the dictionary of tags used by default in the **espadon** package.

Usage

```
dicom.tag.dictionary(add.dict = c("raysearch.tag"))
```

Arguments

add.dict

Vector of the list of additional dictionaries. Put to NULL, if no additional dictionary is requested.

Value

Returns a 3-column dataframe, describing the VR (value representation) and the name of each DICOM TAG.

This dataframe is the fusion of the "nema.tag" dictionary, provided by *nema* [1], with the dictionaries defined in the add.dict vector:

• "raysearch.tag" dictionary is provided by RaySearch laboratories [2]

24 dicom.tag.parser

References

[1] DICOM nema (Online; accessed 2022-02-16). "Current Edition." https://www.dicomstandard.org/current.

[2] Raysearch Laboratories (Online; accessed 2022-04-25). "RAYPLAN 11A, DICOM Conformance Statement." https://www.raysearchlabs.com/4aaf2e/siteassets/raystation-landing-page/dicom-conformance-statements/raystation-pdfs/rsl-d-rs-11a-dcs-en-1.0-2021-05-07-raystation-11a pdf.

Examples

```
str (dicom.tag.dictionary ())
str (dicom.tag.dictionary (NULL))
```

dicom.tag.parser

DICOM TAG parser

Description

the dicom. tag. parser function decodes the content between two DICOM raw data addresses.

Usage

```
dicom.tag.parser(start, stop, VR, endian, dicom.raw.data, try.parse = FALSE)
```

Arguments

Positive integer. Index of the first raw data to parse in the dicom.raw.data.

Positive integer. Index of the last raw data to parse in the dicom.raw.data.

VR Character string, representing the value representation of DICOM data. See DICOM standard.

endian Character string, equal to "little" or "big".

dicom.raw.data Raw vector, representing the binary extraction of the DICOM file.

 $\label{eq:try_parse} \textbf{Boolean. If TRUE, the value, with an undocumented } VR, is considered, as far as$

possible, as a string.

Value

Returns the dicom.raw.data content between the addresses start and stop. Depending on the representation of the value (VR), it can be a character string or a numerical vector.

dicom.to.Rdcm.converter 25

```
dicom.to.Rdcm.converter
```

Conversion of DICOM object into files that can be interpreted by the **espadon** package

Description

The dicom.to.Rdcm.converter function creates, for each DICOM object, a *.Rdcm file usefull for using **espadon** package. Each Rdcm file created is referenced by the date of acquisition of the object (if it is not available, its creation date), the patient's PIN, a reference number, an object number in this reference system, and the object modality (mr, ct, rtstruct...).

Usage

```
dicom.to.Rdcm.converter(
  dcm.files,
  pat.dest.dir,
  update = TRUE,
  tag.dictionary = dicom.tag.dictionary(),
  verbose = TRUE
)
```

Arguments

dcm.files	String vector, representing the list of the full names of the DICOM files of the same patient, or its directory.
pat.dest.dir	Character string, representing the full name of patient directory, which will contain files converted espadon .
update	Boolean. If set to TRUE, and if pat.dest.dir contains previously converted files, these files will keep the same espadon reference frame affectation.
tag.dictionary	Dataframe, by default equal to dicom.tag.dictionary, whose structure it must keep. This dataframe is used to parse DICOM files.
verbose	Boolean. If TRUE, a progress bar indicates the progress of the conversion.

Value

Returns the list of basenames of the created files.

Returns NULL if there are no DICOM files in dcm. files

Note

For each DICOM object, dicom.to.Rdcm.converter creates a *.Rdcm file whose basename is made up of the date of the acquisition (or creation date if previous not found), the patient's PIN, the pseudonym of the frame of reference ("ref1", "ref2"...), the number of the volume object in the directory in this frame of reference ("do1", "do2"...), and the object modality ("mr", "ct", "rtdose", "rtstruct"...).

```
For example: BASE = "20160514_a008e9ac_ref2_do1_mr"
```

26 dicom.viewer

Examples

```
# First, save toy.dicom.raw () raw data to a temporary file for testing.
pat.src.dir <- file.path (tempdir(), "PM_dcm")</pre>
dir.create (pat.src.dir, recursive = TRUE)
dcm.filename <- tempfile (pattern = "PM_rtplan", tmpdir = pat.src.dir,</pre>
                          fileext = ".dcm")
zz <- file (dcm.filename, "wb")</pre>
writeBin (toy.dicom.raw (), zz, size = 1)
close (zz)
\# Create a temporary destination directory where the *.Rdcm file will be saved
pat.dest.dir <- file.path (tempdir(), "PM_Rdcm")</pre>
dicom.to.Rdcm.converter (pat.src.dir, pat.dest.dir, update = TRUE)
dicom.to.Rdcm.converter (dcm.filename, pat.dest.dir, update = TRUE)
list.files (pat.dest.dir)
# Cleaning temporary directories
unlink (pat.src.dir, recursive = TRUE)
unlink (pat.dest.dir, recursive = TRUE)
```

dicom.viewer

DICOM content viewer

Description

the dicom.viewer function displays the data of a DICOM file.

Usage

```
dicom.viewer(
  dcm,
  txt.sep = "\\",
  txt.length = 100,
  tag.dictionary = dicom.tag.dictionary(),
  height = 600,
  width = 900
)
```

Arguments

dcm	String Character representing the list of the full names of a DICOM file or raw data vector of this DICOM file.
txt.sep	String. Used if as.txt = TRUE. Separator of the tag value elements.
txt.length	Positive integer. Used if as.txt = TRUE. Maximum number of letters in the representation of the TAG value.
tag.dictionary	Dataframe, by default equal to dicom.tag.dictionary, whose structure it must keep. This dataframe is used to parse DICOM files.
height, width	Height and width in pixel of the DICOM table.

display.2D.histo

Value

Returns the DICOM file description in a browser window.

See Also

xlsx.from.dcm, dicom.parser

Examples

```
if (interactive ()) dicom.viewer (toy.dicom.raw ())
```

display.2D.histo

Display of a 2D histogram

Description

The display. 2D. histo function displays the density map of a "histo2D" class object.

Usage

```
display.2D.histo(
 histo.2D,
  add = TRUE,
  main = NULL,
  x.lab = NULL,
  y.lab = NULL,
  x.lim = NULL,
  y.lim = NULL,
  bg = "#000000",
  i.rng = NULL,
  display.mode = c("mono.color", "rainbow.color", "line"),
  col = "#ffffff",
  alpha = 1,
  line.pc.levels = c(1, 100),
  line.lwd = 2,
  line.lty = 1
)
```

Arguments

histo.2D	"histo2D" class object.
add	Boolean indicating whether to display the background image.
main	Title of the background image. If $main = NULL$, the title indicates "2D histogram".
x.lab	Label of the x-axis of the background image. If $x.lab = NULL$, this label is $histo.2D$x.file.src$
y.lab	Label of the y-axis of the background image. If $y.lab = NULL$, this label is $histo.2D$y.file.src$.
x.lim	Vector, representing the range of the x-axis.

28 display.2D.histo

y.lim	Vector, representing the range of the y-axis.
bg	Background color of the image. By default, this color is black.
i.rng	Vector of 2 elements giving the minimum and the maximum intensity of the image. If $i.rng = NULL$, then the minimum is 0 and the maximum the maximum density.
display.mode	function display mode. See Details.
col	Color of the color gradient in display.mode = "mono.color", or of the level lines in display.mode = "line". By default, this color is white.
alpha	A number from 0 to 1, indicating the opacity of the image in "rainbow.color" mode. Default alpha = 1 .
line.pc.levels	Vector of level lines in percent of maximum density in display.mode = "line". By default lines 1% and 100% are displayed.
line.lwd	Line thickness of the level lines in display.mode = "line".
line.lty	Type of lines for level lines in display.mode = "line".

Details

The display.mode argument can be set to three values: "mono.color", "rainbow.color", or "line". The 2D histogram graph is displayed by default in "mono.color" mode.

- The "mono.color" mode displays a gradient of the color defined by the col argument, depending on the intensity of \$density.map 2-dimensional array.
- The "rainbow.color" mode makes a display according to the "rainbow" palette, while managing the opacity of the colors.
- The "line" mode draws level lines defined in percent by the line.pc.levels argument.

Value

Returns a display of the density map of histo. 2D. This one must be an object of class "histo2D". See espadon.class for class definitions.

See Also

histo.2D.

display.3D.contour 29

display.3D.contour

Display the 3D contours of the RoI

Description

The display. 3D. contour function performs a 3D display of the selected RoI in the chosen coordinate system.

Usage

```
display.3D.contour(
   struct,
   roi.name = NULL,
   roi.sname = NULL,
   roi.idx = NULL,
   roi.col = NULL,
   roi.print = FALSE,
   roi.lwd = 1,
   roi.cex = 1,
   display.ref = struct$ref.pseudo,
   T.MAT = NULL,
   FoR.axis = FALSE,
   FoR.col = "black"
)
```

Arguments

struct	"struct" class object. See espadon.class for class definitions.
roi.name	Vector of exact names of the RoI in the struct object. By default ${\tt roi.name}$ = NULL. See Details.
roi.sname	Vector of names or parts of names of the RoI in the struct object. By default roi.sname = NULL. See Details.
roi.idx	Vector of indices of the RoI that belong to the struct object. By default roi.idx = NULL. See Details.
roi.col	Color of the RoI. If $roi.col = NULL$ (default), the RoI colors are specified in the $struct$roi.info$.
roi.print	Boolean vector indicating whether to display the pseudonym of the RoI.

30 display.3D.mesh

roi.lwd	Line width of the RoI, by default at 1.
roi.cex	Numeric character expansion factor of RoI name if roi.print = TRUE, defaults to 1.
display.ref	Pseudonym of frame of reference of the display.
T.MAT	"t.mat" class object, created by load.patient.from.Rdcm or load.T.MAT. If T.MAT = NULL, display.ref must be equal to NULL or to struct\$ref.pseudo.
FoR.axis	Boolean or numeric, by default set to FALSE. If FoR.axis = TRUE, the function displays 200 mm length director vectors of the frame of reference. If FoR.axis is numeric, it represent the length in mm of the director vectors.
FoR.col	Color of the frame of reference.

Details

If roi.name, roi.sname, and roi.idx are all NULL, then all of the RoI are selected.

Value

If the concerned regions of interest (RoI) struct exist, it displays the 3D contours of these RoI in the current **RGL** window if it exists, in a new window otherwise.

Examples

display.3D.mesh

3D display of a mesh

Description

The display. 3D. mesh function performs a 3D display of a mesh.

Usage

```
display.3D.mesh(mesh, display.ref = mesh$ref.pseudo, T.MAT = NULL, ...)
```

Arguments

mesh	"mesh" class object, created by the mesh.from.bin function. See espadon.class for class definitions.
display.ref	Character string. Pseudonym of the frame of reference used for display.
T.MAT	"t.mat" class object, created by load.patient.from.Rdcm or load.T.MAT. If T.MAT is NULL, mesh must be displayed in display.ref = mesh\$ref.pseudo.
	Arguments passed to shade3d as color, specular, alpha

display.3D.sections 31

Value

Returns a display of mesh in the current RGL window if it exists, in a new window otherwise.

See Also

mesh.from.bin.

Examples

display.3D.sections Display 3D sections of a patient

Description

The display. 3D. sections function displays transverse, sagittal and frontal views at a point in 3D.

Usage

```
display.3D.sections(
  vol,
  cross.pt = c(0, 0, 0),
  display.ref = vol$ref.pseudo,
  T.MAT = NULL,
  col = grey.colors(10, start = 0, end = 1, alpha = c(rep(0, 1), rep(1, 9))),
  breaks = NULL,
  trans = TRUE,
  sagi = TRUE,
  front = TRUE,
  border = TRUE,
  border.col = "#379DA2"
)
```

32 display.3D.stack

Arguments

vol	"volume" class object to display. See espadon.class for class definitions.
cross.pt	Vector of x, y, z coordinates, representing the cross point of the 3 planes.
display.ref	Character string. Pseudonym of the frame of reference used for display.
T.MAT	"t.mat" class object, created by load.patient.from.Rdcm or load.T.MAT. If T.MAT is NULL, vol must be displayed in display.ref = vol $ref.pseudo$.
col	Vector, representing the color palette of the image. Transparent colors are not represented.
breaks	One of:
	NULL : the minimum and the maximum value of the vol define the range.Vector giving the breakpoints of each color.
trans	Boolean. If TRUE (default), the transverse view is displayed.
sagi	Boolean. If TRUE (default), the sagittal view is displayed.
front	Boolean. If TRUE (default), the frontal view is displayed.
border	Boolean. If TRUE (default), the borders of the planes are displayed
border.col	Color of planes borders

Value

Returns a display of transverse, sagittal and frontal views of vol at cross.pt in the current **RGL** window if it exists, in a new window otherwise. Palette colors are managed by col and breaks.

Examples

display.3D.stack

Display in 3D the selected planes of an espadon class volume

Description

The display. 3D. stack function displays in 3D the requested planes of a "volume" class object.

display.3D.stack 33

Usage

```
display.3D.stack(
  vol,
  k.idx = unique(vol$k.idx[seq(1, vol$n.ijk[3], length.out = 10)]),
  display.ref = vol$ref.pseudo,
  T.MAT = NULL,
  col = grey.colors(10, start = 0, end = 1, alpha = c(rep(0, 1), rep(1, 9))),
  breaks = NULL,
  cube = TRUE,
  border = TRUE,
  ktext = TRUE,
  line.col = "#379DA2",
  line.lwd = 1,
  cex = 1
)
```

Arguments

vol	"volume" class object to display.
k.idx	vector of plane numbers to be displayed, to be chosen in vol\$k.idx. By default k.idx is a vector of 10 uniformly distributed planes in the volume.
display.ref	Character string. Pseudonym of the frame of reference used for display.
T.MAT	"t.mat" class object, created by load.patient.from.Rdcm or load.T.MAT. If T.MAT is NULL, vol must be displayed in display.ref = vol\$ref.pseudo.
col	Vector, representing the color palette of the image. Transparent colors are not represented.
breaks	One of:
	NULL: The minimum and the maximum value of the vol define the range.Vector giving the breakpoints of each color.
cube	Boolean. If TRUE (default), the "volume" edges are displayed.
border	Boolean. If TRUE (default), the borders of the planes defined in k.idx are displayed.
ktext	Boolean. If TRUE (default), the selected plane numbers are displayed.
line.col	Color of cube, planes and texts displayed.
line.lwd	Line width of the border and cube, by default at 1.
cex	Numeric character expansion factor of displayed plan numbers.

Value

Returns a display of the k.idx cutting planes of vol, in the current **RGL** window if it exists, in a new window otherwise. The colors of the palettes are managed by col and breaks.

```
# loading of toy-patient objects (decrease dxyz for better result)
step <- 4
patient <- toy.load.patient (modality = "ct", dxyz = rep (step, 3))
# display o 3 planes</pre>
```

34 display.DVH

display.DVH

Display of a DVH

Description

The display .DVH function displays the Dose Volume Histogram of a "dvh" class object. Y-units are cm^3 .

Usage

```
display.DVH(
  dvh,
  add = FALSE,
  xgrid = TRUE,
  ygrid = TRUE,
  MC.plot = FALSE,
  MC.col = grey.colors(4, rev = TRUE),
  ...
)
```

Arguments

dvh	"dvh" class object.
add	Boolean indicating whether to display the background image.
xgrid	Boolean indicating the display of the x grid.
ygrid	Boolean indicating the display of the y grid.
MC.plot	Boolean. If MC.plot = TRUE, then display. DVH displays, if they exist, the quantile zones (Prob = 0 , .025, .25, .5, .75, .975, 1) of MC DVH variations.
MC.col	Character string, a valid palette with 4 colours corresponding to 100% , 95% , 50% and median of MC data.
	Arguments xlab, ylab, xlim, ylim, main, type, col, lwd, lty and log managed by the plot function.

Value

Returns a plot of the cumulative histogram included in dvh, with its median, and the quantile areas (0%-100%), (2.5%-97.5%) and (25%-75%) of the dvh\$vol variations, if they exist.

See Also

```
display.DVH.pc
```

display.DVH.pc 35

Examples

```
# DVH without MCMC
# loading of toy-patient objects (decrease dxyz and increase beam.nb for
# better result)
step <- 5
patient <- toy.load.patient (modality = c("rtdose", "rtstruct"), roi.name = "",</pre>
                              dxyz = rep (step, 3), beam.nb = 3)
H <- histo.from.roi (patient$rtdose[[1]], patient$rtstruct[[1]], roi.name = "ptv",</pre>
                      breaks = seq (0, 60, by = 2))
DVH <- histo.DVH (H)
display.DVH (DVH)
## Not run:
# loading of toy-patient objects
patient <- toy.load.patient (modality = c("rtdose", "rtstruct"),</pre>
                              roi.name = "gizzard",
                              dxyz = c (2, 2, 2), beam.nb = 3)
# Calculation of the histogram
H <- histo.from.roi (patient$rtdose[[1]], patient$rtstruct[[1]],</pre>
                      roi.name = "gizzard",
                     breaks = seq (0, 60, by = 1), MC = 100)
# DVH
DVH <- histo.DVH (H)
display.DVH (DVH, MC.plot = TRUE, ylim = c (0, 40))
## End(Not run)
```

display.DVH.pc

Display of a cumulative DVH in percent of total volume

Description

The display. DVH. pc function displays the Dose Volume Histogram of "dvh" class object. Y-units are percents of total volume.

Usage

```
display.DVH.pc(
  dvh,
  add = FALSE,
  xgrid = TRUE,
  ygrid = TRUE,
  MC.plot = FALSE,
  MC.col = grey.colors(4, rev = TRUE),
  ...
)
```

Arguments

dvh "dvh" class object. See espadon.class for class definitions.

add Boolean indicating whether to display the background image.

 $display.dV_dx$

xgrid	Boolean indicating the display of the x grid.
ygrid	Boolean indicating the display of the y grid.
MC.plot	Boolean. If MC.plot = TRUE, then display. DVH.pc displays, if they exist, the quantile zones (Prob = 0, .025, .25, .5, .75, .975, 1) of MC DVH variations.
MC.col	Character string, a valid palette with 4 colours corresponding to 100%, 95%, 50% and median of MC data.
	Arguments xlab, ylab, xlim, ylim, main, type, col, lwd, lty and log managed by the plot function.

Value

Returns a plot in percent of total volume of the cumulative histogram included in dvh, with its median, and the quantile areas (0%-100%), (2.5%-97.5%) and (25%-75%) of the dvh\$pcv variations, if they exist.

See Also

```
display.DVH
```

Examples

display.dV_dx

Display of the volume density of a histogram

Description

The display.dV_dx function displays the volume density of a "histo" class object. Y-units are $cm^3.Gy^{-1}$.

Usage

```
display.dV_dx(
  histo,
  add = FALSE,
  xgrid = TRUE,
  ygrid = TRUE,
```

display.histo 37

```
MC.plot = FALSE,
MC.col = grey.colors(4, rev = TRUE),
...
)
```

Arguments

histo	"histo" class object. See espadon.class for class definitions.
add	Boolean indicating whether to display the background image.
xgrid	Boolean indicating the display of the x grid.
ygrid	Boolean indicating the display of the y grid.
MC.plot	Boolean. If MC.plot = TRUE, then display.dV_dx displays, if they exist, the quantile zones (Prob = 0 , .025, .25, .5, .75, .975, 1) of variations in volume density.
MC.col	Character string, a valid palette with 4 colours corresponding to 100% , 95% , 50% and median of MC data.
•••	Arguments xlab, ylab, xlim, ylim, main, type, col, lwd, lty and log managed by the plot function.

Value

Returns a plot of the differential histogram included in histo, with its median, and the quantile areas (0%-100%), (2.5%-97.5%) and (25%-75%) of the histo\$dv_dx variations, if they exist.

See Also

display.histo.

Examples

display.histo

Display of the counts of a histogram

Description

The display.histo function displays the counts of "histo" class object.

38 display.histo

Usage

```
display.histo(
  histo,
  add = FALSE,
  xgrid = TRUE,
  ygrid = TRUE,
  MC.plot = FALSE,
  MC.col = grey.colors(4, rev = TRUE),
  ...
)
```

Arguments

histo	"histo" class object.
add	Boolean indicating whether to display the background image.
xgrid	Boolean indicating the display of the x grid.
ygrid	Boolean indicating the display of the y grid.
MC.plot	Boolean. If MC.plot = TRUE, then display.histo displays, if they exist, the quantile zones (Prob = $0, .05, .25, .5, .75, .95, 1$) of variations in counts.
MC.col	Character string, a valid palette with 4 colours corresponding to 100%, 95%, 50% and median of MC data.
	Arguments xlab, ylab, xlim, ylim, main, type, col, lwd, lty and log managed by the plot function.

Value

Returns a plot of the counts included in histo, with its median, and the quantile areas (0%-100%), (2.5%-97.5%) and (25%-75%) of the histo\$counts variations, if they exist.

See Also

```
display.dV_dx.
```

display.kplane 39

display.kplane

Display of a plane of a volume

Description

The display.kplane function displays the requested plane of a "volume" class object. This function is low-level, used for example in the function display.plane with more intuitive arguments.

Usage

```
display.kplane(
  vol,
  k = vol$k.idx[ceiling(length(vol$k.idx)/2)],
  pt00 = c(0, 0),
  dxy = c(1, 1),
  col = grey.colors(255, start = 0, end = 1),
  breaks = NULL,
  sat.transp = FALSE,
  add = FALSE,
  main = NULL,
  abs.lab = "i",
  ord.lab = "j",
  abs.flip = FALSE,
  ord.flip = FALSE,
  bg = "#000000",
  abs.rng = NULL,
  ord.rng = NULL,
  interpolate = FALSE
)
```

Arguments

vol	"volume" class object to display. See espadon.class for class definitions.
k	Number of the plane to display. By default k is located at mid-plane of the volume.
pt00	Origin point of the displayed plane. By default $pt00 = c (0, 0)$, corresponding to the bottom left of the displayed non-flipped image.
dxy	width and height of a pixel in the plane. If $dxy = c$ (1, 1) (default) abcissa and ordinate correspond to pixel number in the plane.
col	Vector, representing the color palette of the image.
breaks	One of:
	• NULL : the minimum and the maximum value of the vol define the range.
	• Vector giving the breakpoints of each color. Outside values are transparent, leaving the background visible, depending on sat.transp.
sat.transp	Boolean. If TRUE, outside values are transparent, else set to breaks limits colors.
add	Boolean indicating whether to display the background image.
main	Title of the background image. If main = NULL, the title just indicates the value of k.

40 display.kplane

abs.lab	Label of the image abcissa. By default abs.lab = 'i'.
ord.lab	Label of the image ordinate. By default ord.lab = 'j'.
abs.flip	Boolean defaults to FALSE flipping the horizontal axis of the background image.
ord.flip	Boolean defaults to FALSE flipping the vertical axis of the background image.
bg	Background color of the image. By default, this color is black.
abs.rng	Vector of 2 elements indicating the minimum and maximum background image abscissa to display.
ord.rng	Vector of 2 elements indicating the minimum and maximum background image ordinate to display.
interpolate	Boolean, indicating whether to apply linear interpolation to the image.

Value

Returns a display of the k^{th} image plane of vol.

See Also

display.plane.

```
# loading of toy-patient objects (decrease dxyz and increase beam.nb for
# better result)
step <- 5
patient <- toy.load.patient (modality = c("ct","mr", "rtdose"),</pre>
                             dxyz = rep (step, 3), beam.nb = 3)
MR <- patient$mr[[1]]</pre>
CT <- patient$ct[[1]]
D <- patient$rtdose[[1]]</pre>
# display
display.kplane (CT)
display.kplane (MR, k = floor (length(MR$k.idx)*5/8),
                col = grey.colors (256, start = 0, end = 1),
                breaks = seq (0, 500, length.out = 257), bg = "darkblue")
display.kplane (D, k = floor (length(Dk.idx)*3/8),
                col = rainbow (256, s = seq (1, 0, length.out = 256),
                               start = 0, end = 4/6,
                               alpha = seq (0.8, 0, length.out=256),
                               rev = TRUE),
                bg = "darkblue", ord.flip = TRUE, sat.transp = FALSE,
                interpolate = FALSE)
display.kplane (CT, k = floor (length(CT$k.idx)/3), col = pal.RVV (1000),
                breaks = seq(-1000, 1000, length.out = 1001),
                bg = "darkblue", ord.flip = TRUE, interpolate = FALSE)
```

display.legend 41

displa	v.legend	D
u i j D i u	y . ICECIIU	ν

Display of the RoI legend

Description

The display.legend function displays in an image the list of requested RoI and their associated color.

Usage

```
display.legend(
   struct = NULL,
   roi.name = NULL,
   roi.sname = NULL,
   roi.idx = NULL,
   lwd = 1,
   cex = 1,
   displayed.roi.name = NULL,
   bg = "black",
   text.col = "white"
)
```

Arguments

struct	"struct" class object.	
roi.name	Vector of exact names of the RoI in the struct object. By default ${\tt roi.name}$ = NULL. See Details.	
roi.sname	Vector of names or parts of names of the RoI in the struct object. By default roi.sname = NULL. See Details.	
roi.idx	Vector of indices of the RoI that belong to the struct object. By default $\texttt{roi.idx} = \texttt{NULL}$. See Details.	
lwd	Line thickness, defaults to 1	
cex	Font size, default to 1.	
displayed.roi.name		
	Vector. If different from NULL, it represents the replacement names of selected RoI if needed.	
bg	color of the background.	
text.col	color of the legend text.	

Details

roi.name, roi.sname, and roi.idx indicates the RoI to display. If all three are set to NULL, all RoI are selected.

Value

Returns display of the RoI names and their associated color in the active graphics window.

42 display.obj.links

Examples

```
# loading of toy-patient objects
patient <- toy.load.patient (modality = c("rtstruct"), dxyz = c (5, 5, 5))
S <- patient$rtstruct[[1]]
display.legend (struct = S, roi.idx = 2:10, lwd = 2)</pre>
```

display.obj.links

Display patient objects links

Description

The display.obj.links function displays a graph of connections between objects of a patient. The name of the objects corresponds to their modality (ct, mr, rtdose...) followed by their position in their respective lists in the patient list objects. Connected objects are linked by arrows. Objects sharing the same frame of reference have the same color except for objects with warnings, errors or missing planes which are all in grey. Approved objects are circled in green. By default, objects shapes are circles, except rtdose represented as squares.

Usage

```
display.obj.links(
  pat,
  obj.selected = NULL,
  exclusion = NULL,
  square = "rtdose",
  group.by.connected.FoR = TRUE,
  interactive = FALSE,
  random.seed = 314
)
```

Arguments

pat "patient" class object, as loaded using load.patient.from.dicom, load.patient.from.Rdcm

or toy.load.patient.

obj.selected Dataframe (default to NULL) containing the objects already selected, created

by a previous call of display.obj.links for example.

exclusion Vector of patient file modalities that should not be displayed, as for instance

"mr"...

square Vector of patient file modalities that should be enclosed by a square, as for in-

stance c ("ct", "mr")... If NULL no object name is squared.

group.by.connected.FoR

Boolean. If TRUE (default), all objects sharing the same frame of reference or connected by a registration matrix have the same color. If group.by.connected.FoR

=FALSE, only objects sharing the same FoR have the same color.

interactive Boolean. If interactive = TRUE, buttons are available on the graph to get in-

formation about the objects and select or remove them from the data frame of the selected objects. Then simply click on the name of the object on which to apply the chosen action. If interactive = FALSE no interaction possible with

the plot.

display.palette 43

random.seed

Positive Integer or NULL. If random. seed = NULL, the objects are laid out randomly. The layout is otherwise fixed.

Value

The function displays all patient objects, linked by an arrow when they are connected or a line when they belongs to the same DICOM object, and with a color and a shape depending on square, group.by.connected.FoR.

When interactive = TRUE, it returns a dataframe of the selected objects, or NULL if no object is selected.

See Also

get.obj.connection

Examples

```
# loading of toy-patient objects
patient <- toy.load.patient (dxyz = c (5, 5, 5), beam.nb = 1)
display.obj.links (patient, group.by.connected.FoR = FALSE)
display.obj.links (patient, group.by.connected.FoR = TRUE)
display.obj.links (patient, group.by.connected.FoR = TRUE, random.seed=NULL)</pre>
```

display.palette

Display of the color scale of a color palette

Description

The Display.palette function displays the color scale as it is used for representations in espadon functions

Usage

```
display.palette(
  col,
  breaks = NULL,
  factors = NULL,
  override.breaks = FALSE,
  bg = "black",
  new.window = TRUE,
  ylab = ""
)
```

Arguments

col Vector of colors like the ones generated by rainbow, heat.colors, etc.

breaks Vector of breaks for the color palette. It is the usual option for images or dose,

for instance. Its length must be one unit more than col length.

factors Vector containing the labels associated to each col. It should be used for tissue

identification or image segment labelling. Its length must be col length.

44 display.palette

override.breaks

Boolean. When FALSE (by default) ordinates are set to breaks. when TRUE colors are uniformely displayed, and associated breaks set to the correct ordinates for the given colors.

bg Color of the background, seen by transparency for palette having alpha channel.

new.window Boolean. If TRUE, it opens a new window for displaying the palette.

ylab character string. Label of ordinates.

Value

Returns in a new device (if new.window = TRUE), or in the active graphics window (if new.window = FALSE), the palette color defined by col and breaks in priority, or by col and factors.

Note

the breaks are not necessarily evenly spaced. In this case, the colour palette can be represented as the breaks are defined (default option) or by choosing a constant spacing for each colour and displaying the associated abscissa calculated from the breaks (override.breaks = TRUE).

```
## Not run:
# simple example for breaks and factors
display.palette (c ("red", "green", "blue"), breaks = c(0, 1, 3, 7),
                 ylab = "a simple color palette")
display.palette (c ("red", "green", "blue"), breaks = c(0, 1, 3, 7),
                 override.breaks = TRUE)
display.palette (c ("red", "green", "blue"), factors = c ("red", "green", "blue"))
display.palette (c ("gray", "green", "blue"), factors = c (NA, 1, 2))
# for RVV palette, HU range must be [-1000, 1000]
display.palette (pal.RVV (255), breaks = seq (-1000, 1000, length.out=256))
# a palette for dose, for instance
display.palette (rainbow (255, start = 0, end = 4/6, rev = TRUE),
                 breaks = seq (0, 60, length.out = 256))
# black & white palette for CTs or MRs
display.palette (grey.colors (255, start = 0, end = 1),
                 breaks = seq (0, 60, length.out = 256))
# transparency affects colors depending on background (black in first exemple,
# white in the second one)
display.palette (rainbow (255, s = seq (1, 0, length.out = 255),
                          start = 0, end = 4/6,
                          alpha = seq (0.8, 0, length.out = 255), rev = TRUE),
                          breaks = seq (0, 60, length.out=256))
display.palette (rainbow (255, s = seq (1, 0, length.out = 255),
                          start = 0, end = 4/6,
                          alpha = seq (0.8, 0, length.out = 255), rev = TRUE),
                          breaks = seq (0, 60, length.out=256), bg = "white")
## End(Not run)
# colors contracted range using non uniform breaks in the plot window
```

display.plane 45

display.plane

Display the transverse frontal or sagittal view in the patient reference system

Description

The display.plane function displays an overlay of images and RoI closed planar contours on a plane defined by the equations x = constant (sagittal view), or y = constant (frontal view) or z = constant (transverse view) in a frame of reference chosen by the user.

Usage

```
display.plane(
 bottom = NULL,
  top = NULL,
  struct = NULL,
  roi.name = NULL,
  roi.sname = NULL,
  roi.idx = NULL,
  struct.dxyz = c(1, 1, 1),
 display.ref = NULL,
 T.MAT = NULL,
  interpolate = TRUE,
 view.type = c("trans", "front", "sagi"),
 view.coord = 0,
 bg = "#000000",
  abs.rng = NULL,
 ord.rng = NULL,
 bottom.col = grey.colors(255, start = 0, end = 1),
 top.col = rainbow(255, s = seq(1, 0, length.out = 255), start = 0, end = 4/6, alpha =
    seq(0.8, 0, length.out = 255), rev = TRUE),
 bottom.breaks = NULL,
  top.breaks = NULL,
  sat.transp = FALSE,
  struct.lwd = 2,
 main = NULL,
  legend.plot = TRUE,
  legend.shift = 0
)
```

46 display.plane

Arguments

bottom "volume" class object, displayed using bottom.col palette. If bottom = NULL, no bottom image is displayed. "volume" class object, displayed as an overlay, using top.col palette. If top = top NULL, no overlay image is displayed. struct "struct" class object. If NULL, no RoI is displayed. Only RoI of closed planar or point type are displayed. roi.name Vector of exact names of the RoI in the struct object. By default roi.name = NULL. See Details. Vector of names or parts of names of the RoI in the struct object. By default roi.sname roi.sname = NULL. See Details. roi.idx Vector of indices of the RoI that belong to the struct object. By default roi.idx = NULL. See Details. Vector of 3 numbers. Used in case of bottom and top are set to NULL. It represtruct.dxyz sents the virtual steps of a temporary volume created in the display.ref frame of reference, initialized at 1 mm in the 3 directions x, y and z. display.ref Character string. Pseudonym of the frame of reference used for display. If NULL (default), the bottom image FoR, or top image FoR (when no bottom image), or struct FoR (when no volume displayed). T.MAT "t.mat" class object, created by load.patient.from.Rdcm or load.T.MAT. If T.MAT is NULL, bottom, top and struct must have the same frame of reference. interpolate Boolean, indicating whether to apply trilinear interpolation to the bottom and top volumes. If interpolate = FALSE, the values of the nearest voxels are used. When TRUE (by delfault), trilinear interpolation is used. Character string, defining the view to display. It must be set to view.type • "trans" for a transverse view, • "front" for a frontal view or, • "sagi" for a sagittal view. view.coord Numeric vector of the coordinates along the normal vector of the selected view.

Background color of the image. By default, this color is black. bg

Vector of 2 elements indicating the minimum and maximum abscissa to display abs.rng

on the background image.

ord.rng Vector of 2 elements indicating the minimum and maximum ordinate to display

on the background image. bottom.col, top.col

Vectors, representing the palette color of bottom and top.

bottom.breaks, top.breaks

One of:

- NULL : the minimum and the maximum value of bottom or top define the
- Vector giving the breakpoints of each color. Outside values are transparent, leaving the background visible, depending on sat.transp.

When breaks are specified, the number of breaks must be one unit more then the number of colors.

Boolean. If TRUE, outside values are transparent, else set to bottom.breaks or sat.transp top.breaks limits.

display.plane 47

struct.lwd	Line thickness of the RoI contours.
main	Character string. When main different from NULL, it replaces the title, and removes the subtitle and the maximum dose indication if top is of modality rtdose.
legend.plot	Boolean, that indicates whether the RoI legend should be displayed on the image. It is displayed by default.
legend.shift	Numeric. It shifts (in mm) the display of the RoI legend on x-axis.

Details

If roi.name, roi.sname, and roi.idx are all set to NULL, all closed planar or point RoI are selected. If a RoI is not present in the requested plane, the RoI legend won't mention it.

Value

Returns a display of the transverse, sagittal or frontal plane. This plane has the coordinate z = view.coord (transverse), y = view.coord (sagittal) pr x = view.coord (frontal). The display is an overlay of:

- a background image of uniform color bg
- the bottom image if it exists
- the top image if it exists
- the contours of the regions of interest if they exist in the plane considered.

Note

- 1- The main title is given by bottom, the subtitle by top.
- 2- When top is in the "rtdose" modality, the maximum dose is written on the image.

See Also

display.kplane.

```
# loading of toy-patient objects (decrease dxyz and increase beam.nb for
# better result)
step <- 4
patient <- toy.load.patient (modality = c("ct", "mr", "rtstruct", "rtdose"),</pre>
                             roi.name = "",
                              dxyz = rep (step, 3), beam.nb = 3)
CT <- patient$ct[[1]]
MR <- patient$mr[[1]]</pre>
D <- patient$rtdose[[1]]</pre>
S <- patient$rtstruct[[1]]</pre>
display.plane (bottom = CT, top = D, struct = S, view.coord = -30,
               interpolate = FALSE, legend.shift = -80)
# Display of CT in reference frame "ref1" and MR in "ref2"
display.plane (bottom = CT, top = MR, interpolate = FALSE)
# Display of CT and MR in reference frame "ref2"
display.plane (bottom = CT, top = MR, interpolate = FALSE, display.ref ="ref2",
               T.MAT = patient$T.MAT)
```

espadon.class

ESPADON class

Description

ESPADON class

Usage

espadon.class()

Value

Returns a vector of espadon class names.

Note

Each object of a class has specific features that are used to display or process that object.

- the "patient" class includes:
 - \$patient : dataframe providing patient's information as PIN, birth date and gender.
 - \$pat.pseudo: patient's pseudonym, initialized to the patient's PIN of \$patient dataframe.
 - \$description: dataframe describing the patient's DICOM objects: their modality (rtstruct, ct, mr, rtplan ...), the base name of the relevant source file in the patient's directory, the pseudonym of their frame of reference (ref1, ref2 ...), their number of sub-objects, their description if any, their numbers of slices/RoI for all sub-objects, their maximum voxels (for volume sub-objects), and finally the aliases of the sub-objects.
 - \$description.by.reg: list of DICOM objects descriptions that are linked by a transfer matrix.
 - \$T.MAT: list of class "t.mat" containing the information of the transfer matrices to move from one frame of reference to another. See load.T.MAT.
 - \$ct: list of CT, if any. They are named by their \$object.alias See load.obj.from.Rdcm.
 - \$mr: list of MRI, if any. They are formatted like the \$ct.
 - \$rtdose: list of dose matrices. They are formatted like the \$ct.
 - \$rtstruct: list of struct objects.
 - ...any DICOM objects other than the reg files, and those previously mentioned, or any modalities created by **espadon**.
 - \$dicom. dvh: if any, list of DVH computed in rt-dose DICOM files.
- the "t.mat" class includes:
 - \$ref.info: dataframe giving the correspondence between the frame of reference (column \$ref) of the DICOM object (TAG (0020,0052)) and its pseudonym (column \$ref_pseudo).
 - \$reg.info: list of dataframes: the first one gives the PID, birthday, and sex of the patient, the second one gives the name of the source file of transfer matrices.
 - \$matrix.description: dataframe giving the transfer matrix names (column \$t), its source frame of reference (column \$src), the destination frame of reference (column \$dest), and its type (\$type). Note that only the RIGID type is supported.

• \$matrix.list: list of 4X4 transfer matrices. This list contains at least as many Identity matrices as there are ref.pseudo.

A **espadon** object of class "dvh", "fan", "histo", "histo2D", "mesh", "reg", "struct", "t.mat", "undef", "volume" is a list containing at least:

- \$patient: patient's PIN.
- \$patient.bd: patient's birthday.
- \$patient.sex: patient's sex
- \$file.basename: vector of .Rdcm or .dcm file basenames of the object, if it exists
- \$file.dirname: directory including the .Rdcm or .dcm file, if it exists
- \$object.name: name of the object.
- \$object.alias: alias of the object.
- \$frame.of.reference: value of TAG (0020,0052).
- \$ref.pseudo: pseudonym of the \$frame.of.reference
- \$modality: modality of the object (e.g. ct, mr, bin, rtplan..)
- \$description: description of the object.
- \$creation.date: creation date of the object.

If the object was generated from a DICOM file, the list also contains:

- \$object.info: Information of the object. It includes:
 - the SOP ID (value of TAG (0008,0016)),
 - the transfer syntax UID (value of TAG (0002,0010)),
 - the SOP implementation ID (value of TAG (0002,0012)),
 - the SOP type (value of TAG (0008,0008)),
 - the study ID (value of TAG (0020,0010)),
 - the study UID (value of TAG (0020,000D)),
 - the serie UID (value of TAG (0020,000E)),
 - the scanning sequence (value of TAG (0018,0020)),
 - the list of SOP labels (values of TAG (0008,0018)),
 - the dicom source files,
 - the encoding of the content of text tags (values of TAG (0008,0005)) and
 - the number of sub-objects.

if the object is linked to another DICOM object, the list also contains:

- \$ref.object.alias: Alias of the reference object.
- \$ref.object.info: Information of the reference object (not available for mr and ct). It includes:
 - the SOP ID of the reference object,
 - the list of SOP names of the reference object.
- the "dvh" class also includes:
 - \$nb.MC: set to histo\$nb.MC.
 - \$breaks: vector breakpoints.
 - \$mids: vector of cell centers.

• \$mids.unit: Character string, representing the unit of the abcissa of the DVH. For instance, "Gy".

- \$vol: cumulative volume receiving at least the doses defined by \$mids.
- \$pcv: percentage of the total volume receiving at least the doses defined by \$mids.
- if \$nb.MC is different from 0, the arrays MC.vol, MC.pcv and MC.dxyz are added. See histo.DVH.
- the "fan" class also includes:
 - \$origin: the xyz-coordinates of the source point.
 - \$direction: the xyz-coordinates of the main direction of the fan.
 - \$orientation: the xyz-coordinates of the two unit vectors of the plane orthogonal to the \$direction.
 - \$xyz: the xyz-coordinates of the unit vectors of the fan rays
 - \$local: depending on the generation of the fan rays, it can be the spherical coordinates, the deflection angles, the voxel coordinates...
- the "histo" class also includes:
 - \$nb.MC: number of Monte-Carlo simulations
 - \$breaks: vector breakpoints
 - \$mids: vector of cell centers.
 - \$mids.unit: Character string, representing the unit of the abcissa of the histogram. For instance, "Gy".
 - counts: count of voxels whose value is included in the limits defined by \$breaks.
 - dV_dx : differential histogram, expressed in cm^3 by voxel units, at each \$mids.
 - if \$nb.MC is different from 0, the arrays MC.counts, MC.dV_dx and MC.dxyz are added. See histo.from.roi.
- the "histo2D" class also includes:
 - \$nb.pixels: number of elements in the density.map.
 - \$x.file.src: x label. See histo.2D.
 - \$y.file.src: y label. See histo.2D.
 - x.breaks: vector of x-axis breakpoints.
 - y.breaks: vector of y-axis breakpoints.
 - x.mids: vector of x-axis cell centers.
 - y.mids: vector of y-axis cell centers.
 - density.map: array of densities.
 - total.counts: number of counted voxels.
- the "mesh" class also includes:
 - \$nb.faces: set to the number of faces of the mesh.
 - \$mesh: list of 3 elements defining the mesh :\$vb, \$it and \$normals. See mesh.from.bin.
- the "reg" class also includes:
 - \$nb.of.ref: number of transfer matrices.

• \$ref.data: list including the lists of information on transfer matrices, namely: the source frame of reference (\$src), the matrix type (\$type, for example 'RIGID') and the transfer matrix (\$matrix).

- the "rtplan" class also includes:
 - \$approval.status: value of TAG (300E,0002).
 - \$number: sub-object number.
 - \$plan.info: dataframe containing, if they exist,
 - \$label the label for the treatment plan,
 - \$plan. name the name for the treatment plan,
 - \$plan. description description of treatment plan,
 - \$tt.protocol the treatment protocol,
 - \$plan. intent the intent of this plan,
 - \$tt. site describing the anatomical treatment site,
 - \$geometry describing whether RT Plan is based on patient or treatment device geometry.
 - \$presc.dose: dataframe containing, if they exist,
 - \$ref.roi.nb value of TAG (3006,0084),
 - \$dose.ref.nb value of TAG (300A,0012),
 - \$dose.ref.id value of TAG (300A,0013),
 - \$struct. type value of TAG (300A,0014),
 - \$description value of TAG (300A,0016),
 - \$pt. coord value of TAG (300A,0018),
 - \$nominal.prior.dose value of TAG (300A,001A),
 - \$dose. type value of TAG (300A,0020),
 - \$constraint.weight value of TAG (300A,0021),
 - \$deliv.warn.dose value of TAG (300A,0022),
 - \$deliv.max.dose value of TAG (300A,0023),
 - \$targ.min.dose value of TAG (300A,0025), - \$targ.presc.dose value of TAG (300A,0026),
 - \$targ.max.dose value of TAG (300A,0027),
 - \$targ.underdose.vol.frac value of TAG (300A,0028),
 - \$org.risk.full.vol.dose value of TAG (300A,002A),
 - \$org.risk.lim.dose value of TAG (300A,002B),
 - \$org.risk.max.dose value of TAG (300A,002C),
 - \$org.risk.overdose.vol.frac value of TAG (300A,002D)
 - \$fraction.info: dataframe containing, if they exist,
- \$fraction.id the id of the fraction group,
- \$description its description,
- \$planned.frac.nb the total number of treatments (Fractions) prescribed for current fraction group,
- \$frac.pattern.digit.per.day.nb the number of digits in \$frac.pattern used to represent one day,
- \$repeat.frac.cycle.le the number of weeks needed to describe treatment pattern,
- \$frac.pattern the value of TAG (300A,007B) describing treatment pattern every day,
- \$nb.of.beam the number of beams in current fraction group,
- \$beam.dose.meaning the value of TAG (300A,008B) indicating the meaning of Beam Dose,
- \$nb.of.brachy.app the number of brachy application setups in current fraction group.

- \$fraction.beam (in case of beam treatment): dataframe containing, if they exist,
 - \$fraction.id,
 - \$nb.of.frac.planned,
 - \$beam. dose the value of TAG (00A,0084),
 - \$beam. specif.pt the value of TAG (300A,0082),
 - \$beam.meterset the value of TAG (300A,0086),
 - \$beam. type the value of TAG (300A,0090,
 - \$alt.dose the value of TAG (300A,0091),
 - \$alt.type the value of TAG (300A,0092,
 - \$duration.lim the value of TAG (300A,00C5),
 - \$beam. nb the value of TAG (300C,0006) or (300A,00C0),
- \$beam. info (in case of beam treatment): dataframe containing, if they exist,
 - \$beam. nb the value of TAG (300C,0006) or (300A,00C0),
- \$beam.name the value of TAG (300A,00C2),
- \$beam.description the value of TAG (300A,00C3),
- \$beam. type the value of TAG (300A,00C4),
- \$radiation. type the value of TAG (300A,00C6),
- \$high.dose.technique.type the value of TAG (300A,00C7),
- \$treatment.machine.name the value of TAG (300A,00B2),
- \$device.serial.nb the value of TAG (0018,1000),
- \$primary.dosimeter.unit the value of TAG (300A,00B3),
- \$referenced.tolerance.table.nb the value of TAG (300C,00A0),
- \$src.axis.dist the value of TAG (300A,00B4),
- \$referenced.patient.setup.nb the value of TAG (300C,006A),
- \$treatment.delivery.type the value of TAG (300A,00CE),
- \$wedges.nb the value of TAG (300A,00D0),
- \$compensators.nb the value of TAG (300A,00E0),
- \$total.compensator.tray.factor the value of TAG (300A,00E2),
- \$boli.nb the value of TAG (300A,00ED),
- \$blocks.nb the value of TAG (300A,00F0),
- \$total.block.tray.factor the value of TAG (300A,00F2),
- \$final.cumul.meterset.weight the value of TAG (300A,010E),
- \$ctl.pts.nb the value of TAG (300A,0110),
- \$radiation.mass.nb the value of TAG (300A,0302),
- \$radiation.atomic.nb the value of TAG (300A,0304),
- \$radiation.charge.state the value of TAG (300A,0306),
- \$scan.mode the value of TAG (300A,0308),
- \$modulated.scan.mode.type the value of TAG (300A,0309),
- \$virtual.src.axis.dist the value of TAG (300A,030A),
- \$total.wedge.tray.water.equ.thickness the value of TAG (300A,00D7),
- $\t total.compensator.tray.water.equ.thickness the value of TAG (300A,02E3),$
- \$total.block.tray.water.equ.thickness the value of TAG (300A,00F3),
- \$range.shifters.nb the value of TAG (300A,0312),
- \$lateral.spreading.devices.nb the value of TAG (300A,0330),
- \$range.modulators.nb the value of TAG (300A,0340),
- \$fixation.light.azimuthal.angle the value of TAG (300A,0356),
- \$fixation.light.polar.angle the value of TAG (300A,0358).

- \$beam.ctl.pt (in case of beam treatment): list containing, for each beam,
- \$info a data.frame of control points information from DICOM
- \$beam. source the coordinates of the source in the patient frame of reference
- \$beam. direction the coordinates of the beam direction in the patient frame of reference
- \$beam.direction the coordinates of the beam orientation in the patient frame of reference
- \$beam. isocenter the coordinates of the isocenter in the patient frame of reference
- \$spot.map, if they exist, the coordinates of the spots in the patient frame of reference

For the moment, only the rotations of the gantry and the patient support, and the position of the isocenter are taken into account in the calculation of these coordinates.

- \$brachy.info (in case of brachy treatment): dataframe containing, if they exist,
 - \$fraction.id
 - \$nb.of.frac.planned,
 - \$brachy.dose the value of TAG (300A,00A4),
 - \$brachy.nb the value of TAG (300C,000C),
 - \$brachy.specif.pt the value of TAG (300A,00A).
- the "struct" class also includes:
 - \$nb.of.roi: number of regions of interest (RoI).
 - \$thickness: thickness between two consecutive planes of a contour.
 - \$ref.from.contour: reference frame change matrix, from the contour reference frame to the ref.pseudo reference frame
 - \$roi.info: dataframe. Information on RoI contours. It includes the followings columns:
 - -\$number: value of TAG (3006,0084) for the concerned RoI.
 - -\$name: value of TAG (3006,0026) for the concerned RoI.
 - -\$description: value of TAG (3006,0028) for the concerned RoI.
 - -\$generation.algorithm: value of TAG (3006,0036) for the concerned RoI.
 - -\$color: value of TAG (3006,002A) for the concerned RoI.
 - -\$dz: z step between planes for the concerned RoI.
 - -\$roi.pseudo: pseudonym of the RoI \$name. It can be changed by the user.
 - -\$min.x: minimum value x in mm of the RoI. Absent when data = FALSE.
 - -max.x: maximum value x in mm of the RoI. Absent when data = FALSE.
 - -\$min.y: minimum value y in mm of the RoI. Absent when data = FALSE.
 - -\$max.y: maximum value y in mm of the RoI. Absent when data = FALSE.
 - -min.z: minimum value z in mm of the RoI. Absent when data = FALSE.
 - -\$max.z: maximum value z in mm of the RoI. Absent when data = FALSE.
 - -\$vol: volume in mm³ of the RoI. Absent when data = FALSE.
 - -\$Gx: position x in mm of the RoI center of gravity. Absent when data = FALSE.
 - -\$Gy: position y in mm of the RoI center of gravity. Absent when data = FALSE.
 - -\$Gz: position z in mm of the RoI center of gravity. Absent when data = FALSE.
 - -\$continue: boolean, indicating whether the contours are on adjacent planes.
 - \$roi.obs: dataframe. RT RoI observations (sequence TAG (3006,0080)). It includes the followings columns:
 - -\$nb: value of TAG (3006,0082) for the concerned RoI.

- -\$roi.nb: value of TAG (3006,0084) for the concerned RoI.
- -\$label: value of TAG (3006,0085) for the concerned RoI.
- -\$code.value: value of TAG (0008,0100) in the Identification code sequence.
- -\$code.scheme: value of TAG (0008,0102) in the Identification code sequence.
- -\$code.scheme.v: value of TAG (0008,0103) in the Identification code sequence.
- -\$code.meaning: value of TAG (0008,0104) in the Identification code sequence.
- -\$type: value of TAG (3006,00A4) for the concerned RoI.
- -\$interpreter:value of TAG (3006,00A6) for the concerned RoI.
- \$roi.data: exists only if the data is loaded. Contains the list of contour coordinates. The RoI of list number i is that of line i of roi.info. Each element of the list is a list giving the contour information for each plane, namely:
 - \$type: value of TAG (3006,0042).
 - \$pt: dataframe of the coordinates of the contour points.

If the contour is closed (i.e.\$type = "CLOSED_PLANAR"),

then the first point is repeated at the end.

- \$level: contour inclusion level. If this number is even,

the inside of the closed contour belongs to the RoI.

Otherwise, if odd, the inside of the closed contour is excluded from the RoI.

- the "undef" class: is used for DICOM objects that will not be processed further by **espadon** functions. It can also include what the user wants.
- the "volume" class also includes:
 - \$number: sub-object number.
 - \$n.ijk: vector defining the number of indices i, j, k. The product prod(...\$n.ijk) represents the number of voxels in the 3D volume.
 - \$slice. thickness: thickness in mm of a plane.
 - \$min.pixel: minimum value of voxels in the volume.
 - \$max pixel: maximum value of voxels in the volume.
 - \$dxyz: x, y, z steps in mm.
 - \$orientation: value of TAG (0020,0037). Vector, comprising the vectors i and j defining the orientation of the patient with respect to the volume planes.
 - \$xyz0: in the patient frame of reference, position of the first voxel of each plane.
 - \$xyz.from.ijk: transfer matrix of the voxels i, j, k indices to the position x, y, z in mm in the patient's frame of reference.
 - \$k.idx: index of planes in the 3D volume.
 - \$missing.k.idx: Boolean, indicating if k is a continuous sequence of integers.
 - \$cube.idx: 3D volume vertices indices.
 - \$vol3D. data: exists only if the data is loaded. 3D array of the voxel values of the 3D volume.

See Also

toy.load.patient, load.patient.from.dicom, load.patient.from.Rdcm, load.T.MAT histo.DVH, histo.vol, histo.from.roi, histo.from.bin, histo.2D, mesh.from.bin, load.obj.from.Rdcm

```
cat ("espadon class names are:", paste (espadon.class(), collapse = ", "))
```

fan.beam 55

fan.beam

Creation of pyramid fan object with constant angle step.

Description

The fan. beam function creates a "fan" class object containing, among others, the coordinates of the unit director vectors of the rays of rectangular pyramid fan. Rays are uniformly distributed by angle.

Usage

```
fan.beam(
  alpha,
  dalpha,
  orientation = c(0, 0, 1, 1, 0, 0),
  origin = c(0, 0, 0),
  ref.pseudo = "",
  frame.of.reference = "",
  alias = "",
  description = "beam fan"
)
```

Arguments

alpha Positive number specifying the half-angle of the conical beam. dalpha Positive number specifying the step of the angle between the rays of the cone beam. Vector orientation of the pyramid base composed by the 2 orthonormal vectors orientation coordinates. Numeric vector, giving the xyz coordinates of the fan origin. By default c (0, origin ref.pseudo Character string, frame of reference pseudonym of the created object. frame.of.reference Character string, frame of reference of the created object. Character string, \$alias of the created object. alias Character string, describing the the created object. description

Value

Returns a "fan" class object (see espadon.class for class definitions) containing, among others,

- \$xyz: a matrix of 3 columns giving the xyz coordinates of the fan rays.
- \$local: a matrix of 2 columns indicating the deflection angle (in rad) in the main directions defined by orientation.

See Also

fan.planar, fan.sphere, fan.to.voxel.

56 fan.planar

Examples

```
fan <- fan.beam (alpha = 30, dalpha = 1)
head (fan$xyz)
library (rgl)
open3d ()
points3d (fan$xyz)</pre>
```

fan.planar

Creation of pyramid fan object passing through pixels of a plane.

Description

The fan.planar function creates a "fan" class object containing, among others, the coordinates of the unit director vectors of the rays of rectangular pyramid fan. Rays are passing through all pixels of a plane, which represent the pyramid basis. It is for instance useful to compute rt-image.

Usage

```
fan.planar(
  vol,
  k = vol$k.idx[ceiling(length(vol$k.idx)/2)],
  origin = c(0, 0, 0),
  alias = "",
  description = "planar fan"
)
```

Arguments

vol "volume" class object.

k Positive number specifying the plane index that the rays of the fan must cross. By default, k is the central plane.

origin Numeric vector, giving the xyz coordinates of the fan origin. By default c (0, 0, 0).

alias Character string, \$alias of the created object.

description Character string, describing the the created object.

Value

Returns a "fan" class object (see espadon.class for class definitions) containing, among others,

- \$xyz : a matrix of 3 columns giving the xyz coordinates of the fan rays.
- \$local.coords: a list of the ijkt DICOM coordinates of the crossed plane, and the transfer matrix to xyz.from.ijk to compute xyz coordinates in \$ref.pseudo.

See Also

fan.sphere, fan.beam, fan.to.voxel

fan.sphere 57

Examples

```
# loading of toy-patient objects (decrease dxyz and increase beam.nb for
# better result)
step <- 5
patient <- toy.load.patient (modality = c("ct"), dxyz = rep (step, 3))
fan <- fan.planar (patient$ct[[1]], origin = patient$ct[[1]]$xyz0[1,])
head (fan$xyz)
library (rgl)
open3d ()
points3d (fan$xyz)</pre>
```

fan.sphere

Creation of spherical fan object.

Description

The fan. sphere function creates a "fan" class object containing, among others, the coordinates of the unit director vectors of the rays of a spherical fan.

Usage

```
fan.sphere(
  angle = 1,
  method = c("regular", "random"),
  origin = c(0, 0, 0),
  ref.pseudo = "",
  frame.of.reference = "",
  alias = "",
  description = "spherical fan"
)
```

Arguments

	angle	Positive number specifying the angle (or mean angle in case of method = "random") between two nearest vectors.
	method	Requested method of fan calculation from among 'regular' and 'random'. By default, method = regular. See details.
	origin	Numeric vector, giving the xyz coordinates of the fan origin. By default c (0, 0, 0).
	ref.pseudo	Character string, frame of reference pseudonym of the created object.
frame.of.reference		
		Character string, frame of reference of the created object.
	alias	Character string, \$alias of the created object.
	description	Character string, describing the the created object.

58 fan.to.voxel

Details

The "regular" and "random" method are explained by Deserno [1].

• If method = "regular", the returned vectors composing \$xyz matrix are regularly equidistributed at the specified angle.

• If method = "random", the returned vectors composing \$xyz matrix are randomly equidistributed at the specified angle.

Value

Returns a "fan" class object (see espadon.class for class definitions) containing, among others,

- \$xyz: a matrix of 3 columns giving the xyz coordinates of the fan rays.
- \$local: a matrix of 2 columns indicating the polar angle theta (rad) and the azimuthal angle phi (rad) of each ray are added.

References

[1] Deserno, Markus (Online; accessed 2022-08-24). "How to generate equidistributed points on the surface of a sphere." https://www.cmu.edu/biolphys/deserno/pdf/sphere_equi.pdf.

See Also

fan.beam, fan.planar, fan.to.voxel

Examples

```
regular.fan <- fan.sphere (angle = 10)
head (regular.fan$xyz)
random.fan <- fan.sphere (angle = 10, method = "random")
head (random.fan$xyz)
library (rgl)
open3d ()
points3d (regular.fan$xyz)
open3d ()
points3d (random.fan$xyz)</pre>
```

fan.to.voxel

Indices of voxels crossed by a fan

Description

The fan.to.voxel function computes the indices of voxels crossed by a fan. It is useful for retrieving voxel values and voxel indices of a volume (dose or ct) along the fan rays.

Usage

```
fan.to.voxel(vol, fan, restrict = FALSE, vol.value = 1)
```

get.extreme.pt 59

Arguments

vol "volume" class object.

fan "fan" class object created by fan.sphere for example.

restrict Boolean. If TRUE, only the voxels with a value equal to vol.value are taken into account.

vol.value Value of the voxels taken into account, in case of restrict = TRUE

Value

Returns a dataframe of 4 columns. Each line gives:

- column "ray.index": the index (i.e. the row number) of the ray concerned in fan\$dxyz,
- column "vol.index": the index of the voxel crossed in vol\$vol.3Ddata,
- column "l.in": the distance between fan source (i.e. fan\$origin) and the first face of the voxel crossed by the ray,
- column "dl": the distance crossed by the ray in the voxel.

If the rays do not cross any voxel, the dataframe has no row.

See Also

fan.beam, fan.planar, fan.sphere.

Examples

get.extreme.pt

Coordinates of the extreme points

Description

The get.extreme.pt function returns the x, y, z coordinates of the 2 extreme voxels of the rectangular parallelepiped, containing the volume vol. These coordinates are given in the ref. pseudo frame of reference.

Usage

```
get.extreme.pt(vol, ref.pseudo = vol$ref.pseudo, T.MAT = NULL)
```

60 get.ijk.from.index

Arguments

vol "volume" class object.

ref.pseudo Pseudonym of the frame of reference in which you want the coordinates.

T.MAT "t.mat" class object, created by load.patient.from.Rdcm or load.T.MAT. If T.MAT

= NULL, ref.pseudo must be equal to vol\$ref.pseudo.

Value

Returns a dataframe of min and max columns, and x, y and z rows, representing the coordinates of the 2 extreme voxels of the rectangular parallelepiped.

Examples

```
# loading of toy-patient objects
patient <- toy.load.patient (modality = "ct", roi.name = "", dxyz = c (5, 5, 5))
CT <- patient$ct[[1]]
# xyz extreme coordinate
get.extreme.pt (CT)</pre>
```

get.ijk.from.index

Conversion of the indices of a point into ijk vector

Description

The get.ijk.from.index function converts the voxel indices of vol\$vol3D.data (for example, obtained with the function which) into a vector or matrix of DICOM indices i, j, k.

Usage

```
get.ijk.from.index(idx, vol)
```

Arguments

idx Index, or matrix of voxel indices of the array vol\$vol3D.data.

vol "volume" class object.

Value

Returns an i, j, k column matrix of the DICOM indices of the points of vol\$vol3D.data.

See Also

```
get.value.from.ijk, display.kplane
```

get.ijk.from.xyz 61

Examples

get.ijk.from.xyz

Indices relating to the coordinates of the points

Description

The get.ijk.from.xyz function calculates the i, j, k DICOM indices of the points given in the patient's reference frame.

Usage

```
get.ijk.from.xyz(xyz = matrix(c(0, 0, 0), ncol = 3), vol, verbose = FALSE)
```

Arguments

vyz Vector of length 3, corresponding to the x, y, z coordinates (in mm) of a point

in the patient's frame of reference, or 3-column matrix of x, y, z coordinates of

several points.

vol "volume" class object.

verbose Boolean, default to FALSE. If verbose = TRUE, then the xyz coordinates are

printed.

Value

Returns a vector or a matrix of the i, j, k DICOM indices of the x, y, z coordinate points in the patient's frame of reference.

Note

The vector or matrix is made up of real numbers. It is up to the user to make the indices as integer.

The indices of the first voxel vol are 0, 0, 0. WARNING: As i,j,k are DICOM indices, they are not directly related to array indices. To get the value of the vol\$vol3D.data, use the function get.value.from.ijk.

62 get.line

Examples

get.line

Image value along an axis

Description

The get.line function calculates the value of the points of a volume vol along an axis in any direction.

Usage

```
get.line(
  vol,
  origin = c(0, 0, 0),
  orientation = c(1, 0, 0),
  grid = seq(-100, 100, 1),
  interpolate = TRUE
)
```

Arguments

vol "volume" class object.

origin Vector of x, y, z coordinates belonging to the line to extract. If interpolate = FALSE, these coordinates are replaced by the coordinates of the voxel closest to origin.

orientation Directing vector of the line in the vol frame of reference. This vector is internally normalized.

grid Vector, representing the curvilinear coordinates on the line to extract.

interpolate Boolean, default to TRUE. If interpolate = TRUE, a trilinear interpolation of the value of the voxels, relative to the values of adjacent voxels, is performed.

Value

Returns a dataframe, composed of the columns \$x, \$y, \$z, representing the coordinates of the points where the values are taken in vol volume, the column \$s representing the curvilinear abcissa, and the column \$value representing values along \$s.

get.obj.connection 63

Examples

```
# loading of toy-patient objects (decrease dxyz and increase beam.nb for
# better result)
step <- 4
patient <- toy.load.patient (modality = "rtdose", roi.name = "",</pre>
                              dxyz = rep (step, 3), beam.nb = 3)
D <- patient$rtdose[[1]]</pre>
# Dose at maximum dose
origin <- get.xyz.from.index (which.max (D$vol3D.data), D)</pre>
display.plane (bottom = D, view.coord = origin[3],
               bg = "#0000ff")
# Dose profile at x = origin[1] and z = origin[3].
1 <- get.line (vol = D, origin = origin,</pre>
               orientation = c(0, 1, 0), interpolate = FALSE)
plot (1$y, 1$value, type = "1")
grid ()
# Dose profile at y = origin[2] and z = origin[3].
1 <- get.line (D, origin = origin,</pre>
               orientation = c (1, 0, 0), interpolate = FALSE)
plot (1$s, 1$value, type = "1")
grid ()
```

get.obj.connection

List of connections between objects

Description

The get.obj.connection function describes with 4 matrices the different connections between the DICOM objects of the patient.

Usage

```
get.obj.connection(pat)
```

Arguments

pat

"patient" class object, as loaded using load.patient.from.dicom, load.patient.from.Rdcm or toy.load.patient.

Value

Returns a list of 4 named matrices:

- the \$adjacency matrix matrix specifies the source objects that generated the destination objects: the column names correspond to the destinations, and the row names to the sources.
- the \$same.object matrix specifies the elements belonging to the same DICOM object.
- the \$components matrix specifies the objects belonging to the same study.
- the \$same.ref matrix specifies the objects that share the same frame of reference, or with frames of reference linked in T.MAT (by a DICOM reg file for instance)

64 get.plane

See Also

```
display.obj.links
```

Examples

```
# loading of toy-patient objects
patient <- toy.load.patient (dxyz = c (5, 5, 5), beam.nb = 1)
get.obj.connection (patient)
display.obj.links (patient)</pre>
```

get.plane

Extracting a plane from a volume

Description

The get.plane function extracts a plane from a "volume" class object.

Usage

```
get.plane(
  vol,
  origin = c(0, 0, 0),
  plane.orientation = c(1, 0, 0, 0, 1, 0),
  alias = "plane.n",
  xgrid = NULL,
  ygrid = NULL,
  interpolate = TRUE
)
```

Arguments

vol "volume" class object.

origin Vector of x, y, z coordinates, representing the origin of the plane to extract. If

interpolate = FALSE, these coordinates are replaced by the coordinates of the

voxel closest to origin.

plane.orientation

Vector orientation of the plane in the vol frame of reference, composed by the 2 vectors coordinates of the orthonormal basis of the plane. First vector is x-axis,

and second one is y-axis.

alias \$object.alias of the created object.

xgrid Vector, representing the grid of the plane abscissa. See note.

ygrid Vector, representing the grid of the plane ordinates. See note. If ygrid = NULL,

the ordinate is the line intercepting the volume and the step is set to the projec-

tion of vol\$dxyz onto the ordinate orientation.

interpolate Boolean, default to TRUE. If interpolate = TRUE, a trilinear interpolation of the

value of the voxels, relative to the values of adjacent voxels, is performed.

get.plane 65

Value

Returns a "volume" class object, containing only a single plane, at k = 0, in the same frame of reference as vol. This returned object has 2 new fields local.xgrid, and local.ygrid, representing the local grids of the abscissa (columns) and ordinate (rows) of the plane.

Returns NULL if plane doesn't exist.

Note

Determination of axes:

- the x-axis has plane.orientation[1:3] as unit vector.
- the y-axis has plane.orientation[4:6] as unit vector.
- If xgrid is not NULL, origin + x.grid * plane.orientation[1:3] are the coordinates of the points on the x axis.
- If ygrid is not NULL, origin + y.grid * plane.orientation[4:6] are the coordinates of the points on the y axis.
- If xgrid or ygrid are NULL, they are determined to represent as closely as possible the initial volume in the required cut.

```
# loading of toy-patient objects (decrease dxyz for better result)
step <- 4
patient <- toy.load.patient (modality = "mr", dxyz = rep (step, 3))</pre>
MR <- patient$mr[[1]]</pre>
# mid-volume point
mid.point <- apply (get.extreme.pt (MR),1,mean)</pre>
plane <- get.plane (MR, origin = mid.point, interpolate = TRUE)</pre>
display.kplane (plane, interpolate = FALSE)
plane <- get.plane (MR, origin = mid.point, xgrid = seq (-50, 50, 1),
                    ygrid = seq (-50, 50, 1), interpolate = TRUE)
display.kplane (plane, interpolate = FALSE)
# 3 points on the inclined plane
pts <- t ((MR$xyz.from.ijk %*% MR$cube.idx) [1:3 , c (1, 2, 7)])
orientation <- orientation.create (A = pts[1,], B = pts[2,], C = pts[3,])</pre>
origin <- apply (pts, 2, mean)</pre>
plane <- get.plane (MR, origin = origin,</pre>
                     plane.orientation = orientation, interpolate = TRUE)
display.kplane (plane, interpolate = FALSE)
orientation <- orientation.create (A = c (0, 0, 0) , B = c (1, 1, 0),
                                    C = c (-1, 1, 0)
plane <- get.plane (MR, origin = origin,</pre>
                    plane.orientation = orientation, interpolate = TRUE)
display.kplane (plane, interpolate = FALSE)
```

66 get.roi.connection

get.rigid.M	Transfer matrix between two frames of reference	
-------------	---	--

Description

The function get.rigid.M provides, from the T.MAT list created by the functions load.patient.from.Rdcm, load.patient.from.dicom or load.T.MAT, the 4x4 transfer matrix from the FoR (frame o reference) pseudonym src.ref to the FoR pseudonym dest.ref.

Usage

```
get.rigid.M(T.MAT, src.ref, dest.ref)
```

Arguments

T.MAT "t.mat" class object, created by the functions load.patient.from.Rdcm, load.patient.from.dicom or load.T.MAT

src.ref Pseudonym of the source frame of reference

dest.ref Pseudonym of the destination frame of reference

Value

Returns the 4x4 transfer matrix dest.ref from src.ref.

Examples

get.roi.connection List of in

List of inter-connections between RoI

Description

The get.roi.connection function describes the interconnections between Regions of Interest (RoI), from an imaging volume of "cluster" modality, created by struct.clustering.

Usage

```
get.roi.connection(vol)
```

Arguments

vol "volume" class object of "cluster" modality, created by struct.clustering

Value

Returns the list of regions of interest (RoI), where each item in the list gives the inter-connections with other RoI.

get.value.from.ijk 67

See Also

struct.clustering

Examples

get.value.from.ijk

Value of the volume at a selection of DICOM indices

Description

The get.value.from.ijk function calculates the value of a "volume" class object at DICOM indices i, j, k, whether they are integers or not.

Usage

```
get.value.from.ijk(ijk, vol, interpolate = TRUE)
```

Arguments

ijk Vector or 3-column matrix of DICOM indices.

vol "volume" class object.

interpolate Boolean, default to TRUE. If interpolate = TRUE, a trilinear interpolation of the

value of the voxels, relative to the values of adjacent voxels, is performed.

Value

Returns a vector of the values of the volume at the requested DICOM indices.

See Also

get.ijk.from.index.

68 get.value.from.xyz

```
Dmax <- max (D$vol3D.data, na.rm = TRUE)
Dmax
idx <- which (D$vol3D.data >= (Dmax -1) & D$vol3D.data <= (Dmax - 0.2))
ijk <- get.ijk.from.index (idx, D)
get.value.from.ijk (ijk, vol = D, interpolate = FALSE)</pre>
```

get.value.from.xyz

Voxel values on a selection of points

Description

The get.value.from.xyz function calculates the voxel values at the x, y, z coordinate points in the chosen frame of reference.

Usage

```
get.value.from.xyz(
  xyz,
  vol,
  xyz.ref.pseudo = NULL,
  T.MAT = NULL,
  interpolate = TRUE,
  verbose = FALSE
)
```

Arguments

xyz Vector of length 3, corresponding to the x, y, z coordinates (in mm) of a point in

xyz.ref.pseudo frame of reference, or 3-column matrix of x, y, z coordinates

of several points.

vol "volume" class object.

xyz.ref.pseudo ref.pseudo in which the xyz coordinate points are given. This ref.pseudo

must exist in the T.MAT list. If ref.pseudo is NULL then the point with coordinates xyz is considered to be in the patient frame of reference vol\$ref.pseudo.

T.MAT "t.mat" class object, created by load.patient.from.Rdcm, load.patient.from.dicom

or load.T.MAT. If T.MAT = NULL, xyz.ref.pseudo must be equal to vol\$ref.pseudo

or NULL.

interpolate Boolean, default to FALSE. If interpolate = TRUE, a trilinear interpolation of

the value of the voxels, relative to the values of adjacent voxels, is performed.

verbose Boolean, default to FALSE. If verbose = TRUE, then the xyz coordinates are

printed.

Value

Returns a vector of the voxel values at the requested coordinates.

See Also

```
get.xyz.from.index
```

get.volume.from.bin 69

Examples

get.volume.from.bin

Volume selected by binary volume

Description

The get.volume.from.bin function calculates the volume in cm^3 of the selection specified by a "volume" class object of "binary" modality.

Usage

```
get.volume.from.bin(bin)
```

Arguments

bin

"volume" class object, of "binary" modality.

Value

Returns the volume of the binary selection, in cm^3 .

See Also

get.volume.from.roi

70 get.volume.from.roi

```
get.volume.from.roi Volume of a region of interest (RoI)
```

Description

The get.volume.from.roi function extracts the volume $incm^3$ of one or more RoI, from the \$roi.info of the "struct" class object.

Usage

```
get.volume.from.roi(struct, roi.name = NULL, roi.sname = NULL, roi.idx = NULL)
```

Arguments

struct	"struct" class object.
roi.name	Vector of exact names of the RoI in the struct object. By default roi.name = NULL. See Details.
roi.sname	Vector of names or parts of names of the RoI in the struct object. By default roi.sname = NULL. See Details.
roi.idx	Vector of indices of the RoI that belong to the struct object. By default roi.idx = NULL. See Details.

Details

If roi.name, roi.sname, and roi.idx are all set to NULL, all RoI are selected.

Value

Returns a vector of the volumes in cm^3 of the requested RoI.

See Also

get.volume.from.bin, select.names

get.xyz.from.index 71

get.xyz.from.index	Conversion of the indices of a point, into xyz coordinate vector in the
	patient's frame of reference

Description

The get.xyz.from.index function converts the indices of a voxel of vol\$vol3D.data (for example, obtained with the function which) into a vector or matrix of x, y, z coordinates in the patient's frame of reference.

Usage

```
get.xyz.from.index(idx, vol)
```

Arguments

idx Index, or matrix of voxel indices in the array vol\$vol3D.data. The first index of the array is 1.vol "volume" class object.

Value

Returns a column-matrix of coordinates in the patient's reference frame, corresponding to the indices idx.

Examples

grid.equal

Comparison of the grids of two volume objects

Description

The grid.equal function checks that two volumes share the same grid, i.e. the same frame of reference, the same origin point, and the same dx, dy, dz steps.

Usage

```
grid.equal(vol1, vol2)
```

72 histo.2D

Arguments

```
vol1, vol2 "volume" class objects
```

Value

Returns TRUE if the 2 volumes share the same grid.

Examples

histo.2D

2D histograms of 2 volumes

Description

The histo.2D function creates a "histo2D" class object, containing the two-dimensional array of histograms of two "volume" class objects that have the same grid.

Usage

```
histo.2D(
   x.vol,
   y.vol,
   x.breaks = NULL,
   y.breaks = NULL,
   include.outer = TRUE,
   alias = "",
   description = ""
)
```

Arguments

x.vol, y.vol "volume" class objects. The 2 volumes must have the grid (i.e. share the same voxels location).

x.breaks, y.breaks

Vectors giving the breakpoints of x and y axes. See Details.

include.outer Boolean. If include.outer = TRUE, the values out the x.breaks and y.breaks

of each volume are counted in the first and the last cell of the histograms. They

are not taken into account otherwise.

alias Character string, \$alias of the created object

description Character string, describing the created object.

histo.2D 73

Details

The arguments x.breaks and y.breaks represent the scales of the x and y axes of 2D-histogram graph. If they are NULL, the histo. 2D function defaults to 256 cells between the values vol\$min.pixel and vol\$max.pixel.

Value

Returns a "histo2D" class object. This is a list including:

- \$patient: set to x.vol\$patient.
- \$patient.bd: set to x.vol\$patient.bd.
- \$patient.sex: set to x.vol\$patient.sex.
- \$file.basename: set to "".
- \$file.dirname: set to "".
- \$object.name: set to "".
- \$object.alias: alias of the histo2D object.
- \$frame.of.reference: set to x.vol\$frame.of.reference.
- \$ref.pseudo: set to x.vol\$ref.pseudo.
- \$modality: set to "histo2D".
- \$description: description of the histo2D object.
- \$creation.date: set to Sys.Date.
- \$nb.pixels: number of elements in the density.map.
- \$x.file.src: set to x.vol\$object.alias
- \$y.file.src: set to y.vol\$object.alias
- x.breaks: vector of x-axis breakpoints.
- y.breaks: vector of y-axis breakpoints.
- x.mids: vector of x-axis cell centers.
- y.mids: vector of y-axis cell centers.
- density.map: array of densities.
- total.counts: number of counted voxels.

See Also

```
display.2D.histo.
```

74 histo.DVH

```
CT.on.roi <- nesting.roi (CT, S, roi.name = "brain", vol.restrict = TRUE, xyz.margin = c (1, 1, 1), alias = CT$description)

MR.on.CT <- vol.regrid (vol = MR, back.vol = CT.on.roi, interpolate = TRUE, T.MAT = T.MAT, alias = CT$description, description = NULL)

# selection of voxels included in the RoI.

roi.bin <- bin.from.roi (vol = CT.on.roi, struct = S, roi.sname = "brain")

MR.select <- vol.from.bin (MR.on.CT, roi.bin, alias = MR$description)

CT.select <- vol.from.bin (CT.on.roi, roi.bin, alias = CT$description)

# 2D histogram

H2D <- histo.2D (MR.select, CT.select, x.breaks = seq (50, 400, 10), y.breaks = seq (50, 400, 10), alias = "H2D MR CT")

str (H2D)
```

histo.DVH

Cumulative Dose Volume Histogram

Description

The histo. DVH function calculates, for each dose, the volume receiving at least this dose.

Usage

```
histo.DVH(histo, alias = "", description = histo$description)
```

Arguments

histo "histo" class object.

alias Character string, \$alias of the created object.

description Character string, describing the the created object. If the description = NULL

(default value), it will be set to histo\$description.

Value

Returns a "dvh" class object. This is a list including:

- \$patient: set to histo\$patient.
- \$patient.bd: set to histo\$patient.bd.
- \$patient.sex: set to histo\$patient.sex.
- \$file.basename: set to "".
- \$file.dirname: set to "".
- \$object.name: set to "".
- \$object.alias: alias of the dvh object..
- \$frame.of.reference: set to histo\$frame.of.reference.
- \$ref.pseudo: set to histo\$ref.pseudo.
- \$modality: set to "dvh".
- \$description: description of the dvh object. By default, set to histo\$description.
- \$creation.date: set to Sys.Date.

histo.from.bin 75

- \$nb.MC: set to histo\$nb.MC.
- \$breaks: vector breakpoints.
- \$mids: vector of cell centers.
- \$mids.unit: Character string, representing the unit of the abcissa of the DVH. For instance, "Gy", when vol is a rtdose.
- \$vol: cumulative volume receiving at least the doses defined by \$mids.
- \$pcv: percentage of the total volume receiving at least the doses defined by \$mids.
- \$MC.vol: cumulative volume associated with histo\$MC.dV_dx, if it exists.
- \$MC.pcv: percentage of the total volume associated with histo\$MC.dV_dx, if it exists.
- \$MC.dxyz: set to histo\$MC.dxyz, if it exists.

See Also

histo.from.roi, histo.from.bin, histo.vol, display.DVH, display.DVH.pc

Examples

histo.from.bin

Histogram according to a binary

Description

The histo.from.bin function computes the voxels histogram of the selection defined by the binary object sel.bin of a "volume" class object.

Usage

```
histo.from.bin(vol, sel.bin, breaks = NULL, alias = "", description = NULL)
```

76 histo.from.bin

Arguments

vol "volume" class object
sel.bin "volume" class object, of binary modality

breaks Vector giving the breakpoints between histogram cells. If breaks = NULL, the chosen breakpoints are those used by the hist function by default. If breaks are specified, outside values of vol\$vol3D.data are not taken into account.

alias Character string, \$alias of the created object
description Character string, describing the the created object. If the description = NULL (default value),it will be set to vol\$description

Value

Returns a "histo" class object. See histo.vol.

See Also

histo.from.roi, histo.vol, display.histo, display.dV_dx

```
# loading of toy-patient objects (decrease dxyz for better result)
patient <- toy.load.patient (modality = c("ct","rtstruct"), roi.name = "",</pre>
                               dxyz = rep (step, 3))
bin.patient <- bin.from.roi (patient$ct[[1]], struct = patient$rtstruct[[1]],</pre>
                               roi.name = "patient")
# ct histogram in patient volume
H <- histo.from.bin (patient$ct[[1]], sel.bin = bin.patient, breaks = NULL,</pre>
                      alias = "patient_hist")
str(H)
## Not run:
# Skin dose histogram
patient <- toy.load.patient (modality = c("rtdose", "rtstruct"), roi.name = "",</pre>
                               dxyz = c (2, 2, 2), beam.nb = 3)
D <- patient$rtdose[[1]]</pre>
S <- patient$rtstruct[[1]]</pre>
\mbox{\#} Creation of the skin contour of 3 mm
bin.patient <- bin.from.roi (D, struct = S, roi.name = "patient",</pre>
                               alias = "patient")
inverse.patient <- bin.inversion (bin.patient, alias = "inv (patient)")</pre>
expansion <- bin.dilation (inverse.patient, radius = 3,
                             alias = "inv (patient) + 3")
contour.3mm <- bin.intersection (bin.patient, expansion,</pre>
                                   alias = "contour 3 mm")
# Dose histogram in this volume
H <- histo.from.bin (D, sel.bin = contour.3mm, breaks = NULL,
                      alias = "Skin dose")
str(H)
## End(Not run)
```

histo.from.roi 77

nisto.from.roi Histogram according to a RoI

Description

The histo.from.roi function calculates the histogram of the volume voxels belonging to a RoI.

Usage

```
histo.from.roi(
  vol,
  struct,
  roi.name = NULL,
  roi.sname = NULL,
  roi.idx = NULL,
  T.MAT = NULL,
  breaks = NULL,
  MC = NULL,
  sd = c(1, 1, 1),
  offset = c(0, 0, 0),
  alias = "",
  description = NULL
)
```

Arguments

vol	"volume" class object
struct	"struct" class object.
roi.name	Exact name of a RoI in struct object. By default roi.name = NULL. See Details.
roi.sname	Name or part of name of a RoI in struct object. By default roi.sname = NULL. See Details.
roi.idx	Value of the index of a RoI that belong to the struct object. By default roi.idx = NULL. See Details.
T.MAT	"t.mat" class object, created by load.patient.from.Rdcm, load.patient.from.dicom or load.T.MAT. If T.MAT = NULL, struct\$ref.pseudo must be equal to vol\$ref.pseudo.
breaks	Vector giving the breakpoints between histogram cells. If breaks = NULL, the chosen breakpoints are those used by the hist function by default. If breaks are specified, outside values of vol\$vol3D.data are not taken into account.
МС	If different from NULL (default value), number of calculations that will be performed, by Monte-Carlo, by randomly moving the chosen RoI over a random distance, generated according to a normal distribution with mean translation defined by offset and standard deviation sd.
sd	Vector representing the standard deviation of distances in the 3 directions x, y and z.
offset	Vector representing the translation of the RoI in the 3 directions x, y and z.
alias	Character string, \$alias of the created object
description	Character string, describing the the created object. If the description = NULL (default value), it will be set to struct\$roi.info\$roi.pseudo[roi.idx]

78 histo.from.roi

Details

roi.name, roi.sname, and roi.idx must select only one RoI.

Value

Returns "histo" class object. This is a list including:

- \$alias: alias of the histo object.
- \$description: description of the histo object.
- \$breaks: vector breakpoints
- \$mids: vector of cell centers.
- \$mids.unit: Character string, representing the unit of the abcissa of the histogram. For instance, "Gy", when vol is a rtdose.
- counts: count of voxels whose value is included in the limits defined by \$breaks.
- dV_dx: differential histogram, expressed in cm3 by voxel units, at each \$mids.
- MC. counts: array of MC rows. Each row i represents the histogram of the voxels contained in the RoI, whose points have been shifted by \$MC.dxyz[i,].
- MC.dV_dx: array of MC rows. Each row i represents the differential histogram of the voxels contained in the RoI, the points of which have been shifted by \$MC.dxyz[i,].
- MC.dxyz: array of MC rows, representing the offset applied to the RoI.

Note

Using Monte-Carlo can be time consuming for large RoI.

If you only want the result just for a translation, use the arguments MC = 1, sd = 0 and offset = desired translation vector.

See Also

histo.vol, histo.from.bin, display.histo, display.dV_dx

histo.vol 79

	Histogram of a volume	histo.vol
--	-----------------------	-----------

Description

The histo.vol function calculates the voxel values histogram of "volume" class object.

Usage

```
histo.vol(vol, breaks = NULL, alias = "", description = NULL)
```

Arguments

vol "volume" class object.

breaks

Vector giving the breakpoints between histogram cells. If breaks = NULL, the chosen breakpoints are those used by the hist function by default. If breaks are

chosen breakpoints are those used by the hist function by default. If breaks are specified, outside values of vol\$vol3D.data are not taken into account.

alias Character string, \$alias of the created object.

description Character string, describing the the created object. If the description = NULL

(default value), it will be set to vol\$description.

Value

Returns a "histo" class object. This is a list including:

- \$patient: set to vol\$patient.
- \$patient.bd: set to vol\$patient.bd.
- \$patient.sex: set to vol\$patient.sex.
- \$file.basename: set to "".
- \$file.dirname: set to "".
- \$object.name: set to "".
- \$object.alias: alias of the histo object.
- \$frame.of.reference: set to vol\$frame.of.reference.
- \$ref.pseudo: set to vol\$ref.pseudo.
- \$modality: set to "histo".
- \$description: description of the histo object.
- \$creation.date: set to Sys.Date.
- nb.MC: set to 0.
- \$breaks: vector breakpoints
- \$mids: vector of cell centers.
- \$mids.unit: Character string, representing the unit of the abcissa of the histogram. For instance, "Gy", when vol is a rtdose.
- counts: count of voxels whose value is included in the limits defined by \$breaks.
- dV_dx : differential histogram, expressed in cm^3 by voxel units, at each \$mids.

80 load.obj.data

See Also

histo.from.roi, histo.from.bin, display.histo, display.dV_dx

Examples

load.obj.data

Load data of an espadon class object

Description

The load.obj.data function loads all the data of an **espadon** object of class 'struct' or 'volume'.

Usage

```
load.obj.data(obj, tag.dictionary = dicom.tag.dictionary())
```

Arguments

```
obj struct or "volume" class object

tag.dictionary Dataframe, by default equal to dicom.tag.dictionary, whose structure it must keep. This dataframe is used to parse DICOM files in case obj was extracted from DICOM files.
```

Value

Returns the **espadon** object with data \$vol3D.data or \$roi.data

See Also

load.obj.from.dicom and load.obj.from.Rdcm

load.obj.from.dicom

```
CT <- load.obj.data (patient$ct[[1]])
str (CT, max.level = 2)</pre>
```

load.obj.from.dicom Loading an **espadon** object from DICOM files or folder

Description

Loading an **espadon** object from DICOM files or folder.

Usage

```
load.obj.from.dicom(
  dcm.files,
  data = TRUE,
  ref.pseudo = "ref1",
  tag.dictionary = dicom.tag.dictionary(),
  verbose = TRUE
)
```

Arguments

dcm.files	String vector, representing the list of the full names of the DICOM files of the same DICOM object, or its directory.
data	Boolean. Only valid for objects usable by the espadon package, namely ct, mr, rtdose, rtstruct, pt If data = TRUE, either the values of the voxels when modality is (ct, mr, rtdose, pt), or the coordinates of the RoI when modality is rtstruct, are loaded into memory.
ref.pseudo	String, \$ref.pseudo (i.e. pseudonym of the frame of reference) to assign to the loaded object.
tag.dictionary	Dataframe, by default equal to dicom.tag.dictionary, whose structure it must keep. This dataframe is used to parse DICOM files.
verbose	Boolean. If TRUE, a progress bar indicates the progress of the conversion.

Value

Returns an **espadon** object of class "dvh", "histo", "histo2D", "mesh", "rtplan", "struct", "undef" or "volume" depending on the object modality. See espadon.class for class definitions.

See Also

load.obj.data and load.obj.from.Rdcm

82 load.obj.from.Rdcm

```
writeBin (toy.dicom.raw (), zz, size = 1)
close (zz)

# loading of rt-plan object
RTplan <- load.obj.from.dicom (dcm.filename)
str (RTplan)
# Cleaning temporary directory
unlink (pat.dir, recursive = TRUE)</pre>
```

load.obj.from.Rdcm

Loading an **espadon** object from *.Rdcm file

Description

The load.obj.from.Rdcm function loads a DICOM object into memory, creating a list containing the information necessary for its subsequent use with the **espadon** package.

Usage

```
load.obj.from.Rdcm(Rdcm.filename, data = TRUE, nb = NULL)
```

Arguments

Rdcm.filename Character string, representing the full name of a *.Rdcm file created by dicom.to.Rdcm.converter.

Boolean. Only works for objects usable by the **espadon** package, namely ct, mr, rtdose, rtstruct, pt... If data = TRUE, either the values of the voxels when modality is (ct, mr, rtdose), or the coordinates of the RoI when modality is rtstruct, are loaded into memory.

Vector of integers, active only if data = TRUE, and only operating on rtstruct. If nb = NULL, all the RoI of rtstruct are loaded into memory. Otherwise only data of the RoI indices defined by the vector nb are loaded.

Value

Returns an **espadon** object of class "dvh", "histo", "histo2D", "mesh", "rtplan", "struct", "undef" or "volume" depending on the object modality. See espadon.class for class definitions.

See Also

load.obj.data and load.obj.from.dicom

load.patient.from.dicom

Loading patient data from DICOM files

Description

The load.patient.from.dicom function is used to load or pre-load in memory all patient objects from DICOM files.

Usage

```
load.patient.from.dicom(
  dcm.files,
  data = FALSE,
  dvh = FALSE,
  tag.dictionary = dicom.tag.dictionary(),
  verbose = TRUE
)
```

Arguments

dcm.files	String vector, representing the list of the full names of the DICOM files of the same patient, or its directories.
data	Boolean. If data = TRUE, the voxels value of the "volume" class objects, or the coordinates of the RoI (region of interest) of the struct class objects, are loaded into memory.
dvh	Boolean. if dvh = TRUE and if they exist, patient DVH are loaded, for convenience. They are not used as is in espadon package.
tag.dictionary	Dataframe, by default equal to dicom.tag.dictionary, whose structure it must keep. This dataframe is used to parse DICOM files.
verbose	Boolean. If TRUE, a progress bar indicates the progress of the conversion.

Value

Returns an **espadon** object of class "patient", describing the information from dcm.files. See **espadon.class** for a description of the "patient" class.

See Also

dicom.to.Rdcm.converter, load.patient.from.Rdcm, load.obj.data, load.obj.from.dicom, load.obj.from.Rdcm and load.T.MAT.

Examples

```
# First, save toy.dicom.raw () raw data to a temporary file pat.dir for testing.
pat.dir <- file.path (tempdir(), "toy_dccm")</pre>
dir.create (pat.dir, recursive = TRUE)
dcm.filename <- tempfile (pattern = "toyrtplan", tmpdir = pat.dir,</pre>
                           fileext = ".dcm")
zz <- file (dcm.filename, "wb")</pre>
writeBin (toy.dicom.raw (), zz, size = 1)
close (zz)
# loading patient. Here the toy patient ha only a unique rt-plan object
patient <- load.patient.from.dicom (pat.dir, data = FALSE)</pre>
str (patient, max = 2)
# description of object
patient$description
# transfer matrices :
patient$T.MAT
# rt-plan object
str (patient$rtplan[[1]])
# Cleaning temporary directory
unlink (pat.dir, recursive = TRUE)
```

load.patient.from.Rdcm

Loading patient data from *.Rdcm files

Description

The load.patient.from.Rdcm function is used to load or pre-load in memory all patient objects converted in *.Rdcm files.

Usage

```
load.patient.from.Rdcm(dirname, data = FALSE, dvh = FALSE)
```

Arguments

dirname Full paths of the directories of a single patient.

data Boolean. If data = TRUE, the voxels value of the "volume" class objects, or the

coordinates of the RoI (region of interest) of the struct class objects, are loaded

into memory.

dvh Boolean. if dvh = TRUE and if they exist, patient DVH are loaded, for conve-

nience. They are not used as is in **espadon** package.

Value

Returns an **espadon** object of class "patient", describing the information contained in dirname. See **espadon.class** for a description of the "patient" class.

load.Rdcm.raw.data 85

See Also

dicom.to.Rdcm.converter, load.patient.from.dicom, load.obj.data, load.obj.from.dicom, load.obj.from.Rdcm and load.T.MAT.

Examples

load.Rdcm.raw.data

Loading a *.Rdcm file

Description

the load.Rdcm.raw.data function loads the content of a *.Rdcm file.

Usage

```
load.Rdcm.raw.data(Rdcm.filename, address = TRUE, data = TRUE)
```

Arguments

Rdcm.filename Character string, representing the full name of a *.Rdcm file created by dicom.to.Rdcm.converter.

address boolean. If TRUE, a dataframe with the address of the tags in the raw DICOM data is returned.

boolean. If TRUE, the DICOM information are returned as an R list.

Value

Returns a list containing the information, converted by **espadon**, of a DICOM object..

See Also

dicom.to.Rdcm.converter, load.obj.from.Rdcm.

86 load.T.MAT

Examples

load.T.MAT

Loading of information about transfer matrices between frames of reference of patient objects.

Description

The load.T.MAT function lists all the frames of reference of the objects included in the patient directory. It concatenates all the information of the reg matrices of a directory, creating, among other things, a list of 4x4 transfer matrices between frames of reference.

Usage

```
load.T.MAT(dirname)
```

Arguments

dirname

Directory containing a patient's DICOM objects.

Value

Returns a "t.mat" class object. It is a list that includes:

- \$ref.info: dataframe giving the correspondence between the frame of reference (column \$ref) of the DICOM object (TAG (0020,0052)) and its pseudonym (column \$ref_pseudo).
- \$reg.info:list of dataframes: the first one gives the PID, birthday, and sex of the patient, the second one gives the name of the source file of transfer matrices.
- \$matrix.description: dataframe giving the transfer matrix names (column \$t), its source frame of reference (column \$src), the destination frame of reference (column \$dest), and its type (\$type). Note: only the RIGID type is supported.
- \$matrix.list: list of 4X4 transfer matrices. This list contains at least as many Identity matrices as there are ref.pseudo.

mesh.from.bin 87

Examples

mesh.from.bin

Creation of a mesh according to a binary volume

Description

The mesh. from. bin function creates a mesh class object from a volume object of "binary" modality.

Usage

Arguments

bin "volume" class object of "binary" modality.

alias Character string, \$alias of the mesh defining the \$alias of the created object.

tol Tolerance in mm, applied for mesh simplification. See vcgClean.

smooth.iteration

Number of iterations applied in the smoothing algorithm. See vcgSmooth.

smooth.type character: select smoothing algorithm. Available are "taubin", "laplace", "HClaplace", "fujiLaplace", "angWeight" (and any sensible abbreviations). By default, set to "taubin". See vcgSmooth.

smooth.lambda numeric: parameter for Taubin smooth. See vcgSmooth.

88 mesh.from.bin

smooth.mu numeric: parameter for Taubin smooth. See vcgSmooth.

smooth.delta numeric: parameter for Scale dependent laplacian smoothing (see reference below).and maximum allowed angle (in radians) for deviation between normals Laplacian (surface preserving). See vcgSmooth.

verbose Boolean, by default set to FALSE. Allows you to inhibit comments.

Value

Returns a "mesh" class object. This is a list including the following 6 elements:

- \$patient: set to bin\$patient.
- \$patient.bd: set to bin\$patient.bd.
- \$patient.sex: set to bin\$patient.sex.
- \$file.basename: set to "".
- \$file.dirname: set to "".
- \$object.name: set to "".
- \$object.alias: set to the alias argument of the function.
- \$frame.of.reference: set to bin\$frame.of.reference.
- \$ref.pseudo: set to bin\$ref.pseudo.
- \$modality: set to "mesh".
- \$description: By default, set to paste (bin\$object.alias, "mesh").
- \$creation.date: set to Sys.Date.
- \$nb. faces: set to the number of faces of the mesh.
- \$mesh: list of 3 elements defining the mesh:
- \$vb: array made up of the generalized coordinates (x, y, z, 1) of the vertices of the triangles.

There are as many columns as there are vertices.

- \$it: array of the 3 indices of the vertices forming a triangle, arranged by column.

There are as many columns as there are triangles in the mesh.

- normals: array made up of the generalized coordinates (x, y, z, 1) of the normal vectors of each triangle. There are as many columns as there are vertices.

Note

To compute the mesh, all NA voxels of the binary volume bin are set to FALSE. If all voxels are equal to FALSE, the function returns the code NULL.

See Also

vcgSmooth

mesh.in.new.ref

```
S <- patient$rtstruct[[1]]
# creation of the patient mesh
bin <- bin.from.roi (CT, struct = S, roi.name = "patient")
mesh.patient <- mesh.from.bin (bin, alias = "patient", verbose = FALSE)
str (mesh.patient)</pre>
```

mesh.in.new.ref

Change of frame of reference of a mesh

Description

The mesh.in.new.ref function allows you to change the frame of reference of a mesh.

Usage

```
mesh.in.new.ref(mesh, new.ref.pseudo, T.MAT, alias = "", description = NULL)
```

Arguments

mesh "volume" class object.

new.ref.pseudo pseudonym of the frame of reference in which the mesh should be located. This new.ref.pseudo must exist in the T.MAT list.

T.MAT "t.mat" class object, created by load.patient.from.Rdcm, load.patient.from.dicom, load.T.MAT or ref.add.

alias Character string, \$alias of the created object.

Character string, describing the created object. If description = NULL (default value), it will be that of the mesh.

Value

Returns "mesh" class object in the new frame of reference new.ref.pseudo.

90 mesh.spheric.proj

mesh.repair

Repair of a mesh

Description

The mesh.repair function repairs holes in a mesh class object.

Usage

```
mesh.repair(mesh, verbose = TRUE)
```

Arguments

mesh

"mesh" class object.

verbose

Boolean, by default set to FALSE. Allows you to inhibit comments.

Value

Returns a mesh, repaired by removing degenerated triangles and filling holes.

Examples

mesh.spheric.proj

Adding spherical coordinates to a mesh

Description

The mesh.spheric.proj function adds latitude and longitude coordinates to a mesh. These features map the mesh surface to a sphere. Latitude and longitude are computed using the heat diffusion approach explained by *Brechbühler and al* [1].

Usage

```
mesh.spheric.proj(mesh, verbose = TRUE)
```

mesh.spheric.proj 91

Arguments

mesh "mesh" class object.

verbose Boolean, by default set to FALSE. Allows you to inhibit comments.

Value

returns a "mesh" class object in which mesh contains Lat and 1on evaluated at vertices. The function allows to have a parameterized surface for later computations as curvature or shape index, hence, nor the surface, nor the angles are preserved. In the DICOM frame of reference, latitude goes along Z axis (from feet = -1 to head = +1) and longitude turns counter clockwise (from -1 to +1).

Note

This funtion is time consuming.

References

[1] Brechbuhler C, Gerig G, Kubler O (1995). "Parametrization of Closed Surfaces for 3-D Shape Description." *Computer Vision and Image Understanding*, **61**(2), 154-170. ISSN 1077-3142, doi:10.1006/cviu.1995.1013.

```
# loading of toy-patient objects (decrease dxyz for better result)
patient <- toy.load.patient (modality = c("ct", "rtstruct"), roi.name = "",</pre>
                               dxyz = rep (step, 3))
CT <- patient$ct[[1]]
S <- patient$rtstruct[[1]]</pre>
#creation of the patient mesh
bin <- bin.from.roi (CT, struct = S, roi.name = "patient")</pre>
m.patient <- mesh.from.bin (bin)</pre>
m.skin <- mesh.repair (m.patient, verbose = FALSE)</pre>
m.proj <- mesh.spheric.proj (m.skin, verbose = FALSE)</pre>
library (rgl)
col <- hcl.colors (12, "Blue-Red 3")</pre>
open3d()
shade3d (m.proj$mesh, meshColors = "vertices",
         color = col[round ((m.proj$mesh$Lat/2 + 0.5) * 11) + 1],
         specular = "#404040")
open3d()
shade3d (m.proj$mesh, meshColors = "vertices",
         color = col[round ((m.proj$mesh$Lon/2 + 0.5) * 11) + 1],
         specular = "#404040")
```

92 nesting.bin

nesting.bin

Restrict volume to a binary selection

Description

The nesting .bin function restricts a "volume" class object to the rectangular parallelepiped circumscribed to the selected voxels.

Usage

```
nesting.bin(
  vol,
  sel.bin,
  alias = "",
  description = NULL,
  xyz.margin = c(0, 0, 0)
)
```

Arguments

"volume" class object, containing data to restrict.

sel.bin "volume" class object, of "binary" modality, specifying the selected voxels.

Character string, \$alias of the created object.

Character string, describing the the created object. If description = NULL, it will be paste (vol\$description, "restricted to", sel.bin\$description).

xyz.margin Vector of length 3, representing the distances in mm to be added to the x, y and z directions of the rectangular parallelepiped circumscribed to the voxels selected in sel.bin, in the cutting planes frame of reference. By default xyz.margin = c (0, 0, 0).

Value

Returns a "volume" class object, in which 3D volume is limited to the rectangular parallelepiped circumscribed to the voxels selected by sel.bin, increased by the requested margins.

nesting.cube 93

nesting.cube	Restriction of a volume to a rectangular parallelepiped	

Description

The nesting. cube function restricts or increases a volume to the rectangular parallelepiped defined by its 2 extreme vertices.

Usage

```
nesting.cube(vol, pt.min, pt.max, alias = "", description = NULL)
```

Arguments

vol	"volume" class object.
pt.min	minimum x, y, z coordinates of the rectangular parallelepiped vertex.
pt.max	maximum x, y, z coordinates of the rectangular parallelepiped vertex.
alias	Character string, \$alias of the created object.
description	Character string, describing the the created object. If the description = NULL (default value), it will be set to vol\$description.

Value

Returns a "volume" class object, in which 3D volume is restricted or increased to be circumscribed to the requested rectangular parallelepiped. If the created volume exceeds the initial volume, new voxels are set to NA.

See Also

add.margin.

94 nesting.roi

οI	
----	--

Description

The nesting.roi function restricts a "volume" class object to the rectangular parallelepiped circumscribed to the chosen RoI.

Usage

```
nesting.roi(
  vol,
  struct,
  roi.name = NULL,
  roi.sname = NULL,
  roi.idx = NULL,
  xyz.margin = c(0, 0, 0),
  vol.restrict = FALSE,
  T.MAT = NULL,
  alias = "",
  description = NULL
)
```

Arguments

vol	"volume" class object
struct	"struct" class object.
roi.name	Vector of exact names of the RoI in the struct object. By default roi.name = NULL. See Details.
roi.sname	Names or parts of names of the RoI in the struct object. By default roi.sname = NULL. See Details.
roi.idx	Index of the RoI that belong to the struct object. By default roi.idx = NULL. See Details.
xyz.margin	Vector of length 3, representing the distances in mm to be added to the x, y and z directions of the rectangular parallelepiped circumscribed to the chosen RoI, in the cutting planes frame of reference. By default xyz.margin = $c(0, 0, 0)$.
vol.restrict	Boolean. If vol.restrict = TRUE, the rectangular parallelepiped circumscribed to the chosen RoI, enlarged by xyz.margin cannot exceed the initial volume.
T.MAT	"t.mat" class object, created by load.patient.from.dicom, load.patient.from.Rdcm or load.T.MAT. If T.MAT = NULL, struct\$ref.pseudo must be equal to vol\$ref.pseudo.
alias	Character string, \$alias of the created object.
description	Character string, describing the the created object. If description = NULL, it will be that of the vol.

Details

If $\verb"roi.name", \verb"roi.sname", and \verb"roi.idx" are all set to NULL, all RoI are selected.$

orientation.create 95

Value

Returns a "volume" class object, in which 3D volume is limited to the rectangular parallelepiped circumscribed to the chosen RoI, increased by the requested margins.

Examples

```
# loading of toy-patient objects (decrease dxyz for better result)
step <- 4
patient <- toy.load.patient (modality = c("ct", "rtstruct"),</pre>
                              roi.name = "brain", dxyz = rep (step, 3))
CT <- patient$ct[[1]]</pre>
S <- patient$rtstruct[[1]]</pre>
CT.brain <- nesting.roi (CT, S, roi.sname = "brain")</pre>
CT.brain.with.margin <- nesting.roi (CT, S, roi.sname = "brain",
                                           xyz.margin = c (10,10,10))
\mbox{\tt\#} display at the center of gravity of the cerebellum \mbox{\tt Gz}
Gz <- S$roi.info [grep("^brain",S$roi.info$roi.pseudo),]$Gz</pre>
display.plane (bottom = CT.brain, view.coord = Gz,
               struct = S, bottom.col = pal.RVV (1000),
               bottom.breaks = seq (-1000, 1000, length.out = 1001),
               bg = "#00ff00", interpolate = FALSE, legend.shift = -20)
display.plane (bottom = CT.brain.with.margin,view.coord = Gz,
               struct = S, bottom.col = pal.RVV (1000),
               bottom.breaks = seq(-1000, 1000, length.out = 1001),
               bg = "#00ff00", interpolate = FALSE, legend.shift = -20)
```

orientation.create

Creation of orientation

Description

The orientation.create function creates the orientation vectors of a plane:

- from 3 points A, B and C (see details),
- or from 2 vectors B and C, resp. defining x and y-axis (see details),
- or from 2 points A, B defining x-axis, and the normal vector to the plane (see details),
- or from a vector B defining x-axis, and the normal vector to the plane (see details).

Usage

```
orientation.create(A = c(0, 0, 0), B = NULL, C = NULL, normal = NULL)
```

Arguments

A	Vector of the x, y and z coordinates of point A, by default equal to $c(0,0,0)$ in the case where B and C are vectors.
В	Vector of x, y and z coordinates of point or vector B.
С	Vector of x, y and z coordinates of point or vector C.
normal	Vector of x, y and z coordinates of normal vector.

96 pal.RVV

Details

When using B and C, B-A define the x-axis unit vector. The unit vector of the y-axis is orthonormal to the x-axis, coplanar with A, B and C, and in the direction of A to C.

When using B and normal, the unit vector of the x-axis is orthonormal to the normal vector, in the direction of A to B. The unit vector of the y-axis is defined so as to constitute a direct orthonormal basis with the unit vector of the x-axis and the normal vector of the plane.

Value

Returns the orientation of the plane. That means the concatenation of 2 vectors, defining an orthonormal basis of the plane.

Examples

```
A <- c (-29.93, 18.85, 4.34)
B <- c (28.73, 15.36, 4.46)
C <- c (1.53, 75.21, 13.51)
orientation.create (A, B, C)
```

pal.RVV

Conversion of Hounsfied Units to Realistic Volume Vizualization colors

Description

The RVV.pal function produces a color palette where Hounsfield Units in the range -1000 HU to 1000 HU are converted into realistically colorized virtual anatomy (for use with CT), developed by *J.C. Silverstein and al* [1]

Usage

```
pal.RVV(n, alpha = NULL, min.col = "#000000", max.col = "#FFFFFF")
```

Arguments

n Integer, number of colors to be in the palette

alpha Vector representing the opacity, in the range of 0 (transparent) to 1 (opaque). If

alpha = NULL (default), all colors are opaque, and no alpha channel is added to

the colors.

min.col, max.col

respectively the color below -1000HU (by default, black, i.e. "#000000") and

above +1000HU (by default, white, i.e. "#FFFFFF")

Value

Returns a vector of colors of size n.

References

[1] Silverstein JC, Parsad NM, Tsirline V (2008). "Automatic perceptual color map generation for realistic volume visualization." *Journal of Biomedical Informatics*, **41**(6), 927-935. ISSN 1532-0464, doi:10.1016/j.jbi.2008.02.008.

Rdcm.inventory 97

Examples

```
pal <- pal.RVV (256)

image (x = seq (-1000, 1000, length.out = 1024), y = 1,
    z = matrix (seq (-1100, 1100, length.out = 1024), ncol = 1),
    col = pal,
    main = "Realistic Volume Vizualization colors")</pre>
```

Rdcm.inventory

Inventory of espadon objects from Rdcm files

Description

The Rdcm.inventory function creates, from Rdcm files in a patient's directory, a dataframe describing objects.

Usage

```
Rdcm.inventory(dirname)
```

Arguments

dirname

Character string, representing the full name of patient directory, including Rdcm files.

Value

Returns a dataframe, providing information of DICOM objects.

98 ref.add

ref.add

Adding a frame of reference in T.MAT

Description

The ref. add function adds the transfer matrices from or to a new frame of reference defined from 2 unit vectors and an origin point.

Usage

```
ref.add(
    src.ref,
    orientation = c(1, 0, 0, 0, 1, 0),
    origin = c(0, 0, 0),
    new.ref.pseudo = "newref",
    T.MAT = NULL
)
```

Arguments

Character string, pseudonym of the frame of reference in which the orientation vector and the origin point origin are defined.

Orientation

Vector of 6 or 9 elements, composed of the coordinates of the 2 orthonormal vectors (i, j), or of the 3 orthonormal vectors (i, j, k) of the new coordinate system, in the src.ref frame of reference.

Origin

Vector of the x, y, z coordinates of the origin point of the new frame of reference in the src.ref frame of reference. Default to c (0, 0, 0).

New.ref.pseudo

Character string, pseudonym of the new frame of reference to add.

T.MAT

"t.mat" class object created by load.patient.from.dicom, load.patient.from.Rdcm or load.T.MAT. If T.MAT = NULL, then only the link between src.ref and new.ref.pseudo is computed.

Value

Returns a "t.mat" class object, which contains the transfer matrices from or to new.ref.pseudo pseudonym of the new frame of reference. If the T.MAT is NULL, then the returned object will contain only 4 matrices: "src.ref<-src.ref", "src.ref<-new.ref.pseudo", "new.ref.pseudo<- new.ref.pseudo<-, "new.ref.pseudo<- ref".

Returns a NULL if orientation is not well defined.

See Also

ref.cutplane.add, ref.remove.

ref.cutplane.add 99

ref.cutplane.add

Adding volume's cutting planes frame of reference in T.MAT

Description

The ref.cutplane.add function adds in T.MAT the transfer matrices from or to volume's cutting planes frame of reference.

Usage

```
ref.cutplane.add(
  vol,
  origin = vol$xyz0[1, ],
  ref.cutplane = paste0(vol$ref.pseudo, "m"),
  T.MAT = NULL
)
```

Arguments

vol "volume" class object.

origin Vector of the x, y, z coordinates of the origin point of the cut planes frame of

reference.

ref.cutplane Name of the volume's cutting planes frame of reference. By default ref.cutplane

= paste0 (vol\$ref.pseudo,".m").

T.MAT "t.mat" class object created by load.patient.from.dicom, load.patient.from.Rdcm

or load.T.MAT. If T.MAT = NULL, then only the link between vol\$ref.pseudo

and ref. name is established.

Value

Returns a "t.mat" class object, which contains the transfer matrices from or to volume's cutting planes frame of reference. If the T.MAT is NULL, then the returned object will contain only 4 matrices: "src.ref<-src.ref", "src.ref<-ref.cutplane", "ref.cutplane", "ref.cut

See Also

ref.add, ref.remove.

```
# loading of toy-patient objects
patient <- toy.load.patient (modality = "mr", roi.name = "", dxyz = c (4, 4, 4))
MR <- patient$mr[[1]]
MR$xyz.from.ijk
# creation of t.mat, containing the transfer matrix to the frame of reference</pre>
```

100 ref.remove

```
# of the MR cutting planes
t.mat <- ref.cutplane.add (MR)

# Change of frame of reference
MR.m <- vol.in.new.ref (MR, paste0 (MR$ref.pseudo, "m"), t.mat)
MR.m$xyz.from.ijk</pre>
```

ref.remove

Deletion of a frame of reference in T.MAT

Description

The ref. remove function removes the management of a frame of reference in T.MAT.

Usage

```
ref.remove(ref.name, T.MAT)
```

Arguments

ref.name Character string, pseudonym of the frame of reference to delete.

T.MAT "t.mat" class object in which the ref. name frame of reference is to be deleted.

Value

Returns a "t.mat" class object, which no longer contains transfer matrices from or to the ref.pseudo ref.name.

See Also

ref.add

ref.srctodest.add 101

ref.srctodest.add	Linking two	existing frames	of reference	in TMAT
rer.srctodest.add	Linking two	existing traines	oi reierence	e in L.WIAT

Description

The ref.srctodest.add function links the source frame of reference with the destination frame of reference.

Usage

```
ref.srctodest.add(src.ref, dest.ref, TM = diag(4), T.MAT = NULL)
```

Arguments

src.ref	Character string, pseudonym of the source frame of reference.
dest.ref	Character string, pseudonym of the destination frame of reference.
TM	4x4 tansfert matrix for moving from src.ref to dest.ref.
T.MAT	"t.mat" class object created by load.patient.from.dicom, load.patient.from.Rdcm or load.T.MAT. If T.MAT = NULL, then only the link between src.ref and dest.ref is established.

Value

Returns a "t.mat" class object, which contains the transfer matrices from or to dest.ref pseudonym of the new frame of reference. If the T.MAT is NULL, then the returned object will contain only 4 matrices: "src.ref<-src.ref", "src.ref<-dest.ref", "dest.ref<- dest.ref", "dest.ref<-src.ref".

See Also

ref.add, ref.cutplane.add, ref.remove.

102 rt.chi.index

rt.chi.index Chi index 2D - 3D

Description

The rt.chi.index function computes the local or global Chi index from a measurement and a reference. These latter are "volume" class objects containing one (2D) or several planes (3D).

Usage

```
rt.chi.index(
  vol,
  vol.ref,
  abs = TRUE,
  vol.max = vol.ref$max.pixel,
  dose.th = 0.02,
  delta.r = 3,
  analysis.th = 0.05,
  local = FALSE,
  local.th = 0.3,
  alias = "",
  description = NULL
)
```

Arguments

vol	"volume" class object, which represents the measured volume.
vol.ref	"volume" class object, which represents the reference volume.
abs	Boolean. If TRUE (default), the absolute value of Chi is computed.
vol.max	Positive number, by default equal to the maximum value of the reference volume. See Details.
dose.th	Positive number, in percent, used to determine the dose difference criterion. See Details.
delta.r	Positive number, in mm. Distance difference criterion.
analysis.th	Positive number, in percent. Only the voxels whose value are greater than or equal analyse.th \star vol.max are processed.
local	Boolean. If local = FALSE (default), a global Chi index is computed, and a local Chi index otherwise.
local.th	Positive number, in percent. Local threshold, only used if local = TRUE. See Details.
alias	Character string, \$object.alias of the created object.
description	Character string, describing the created object. If description = NULL (default value), it will be set to Chi index setup.

rt.chi.index 103

Details

The Chi index of a voxel n was defined by Bakai and al [1]. It is computed from the formulae:

$$\chi_n = \frac{D_i - Dref_n}{\sqrt{\Delta D^2 + \Delta r^2 \cdot \|\nabla Dref_n\|^2}}$$

If abs = TRUE, the used formulae is:

$$\chi_n = \frac{|D_i - Dref_n|}{\sqrt{\Delta D^2 + \Delta r^2 \cdot \|\nabla Dref_n\|^2}}$$

with D_i the measured dose at voxel i, $Dref_n$ the reference dose at voxel n, $\nabla Dref_n$ the gradient of reference dose at voxel n, Δr the distance difference criterion equal to delta.r, and ΔD the distance difference criterion at voxel n defined as follows:

- If local = FALSE a global Chi index is computed and $\Delta D = dose.th \cdot vol.max$.
- If local = TRUE, then $\Delta D = dose.th \cdot Dref_n$ when $Dref_n \geq local.th \cdot vol.max$, and $\Delta D = dose.th \cdot local.th \cdot vol.max$ otherwise.

Value

Returns a "volume" class object (see espadon.class for class definitions). The \$vol3D.data field represents the Chi index. Two fields are added: the \$setup field recalls the calculation setup, and the \$chi.info field details the rate of absolute values of the Chi index below 1, above 1.2 and 1.5, the max and the mean Chi index.

References

[1] Bakai A, Alber A, Nüsslin F (2003). "A revision of the Gamma-evaluation concept for the comparison of dose distributions." *Physics in Medicine and Biologys*, **48**(21), 3543–3553.

See Also

rt.gamma.index

104 rt.gamma.index

rt.gamma.index

Gamma index 2D - 3D

Description

The rt.gamma.index function computes the local or global Gamma index from a measurement and a reference. These latter are "volume" class objects containing one (2D) or several planes (3D).

Usage

```
rt.gamma.index(
  vol,
  vol.ref,
  vol.max = vol.ref$max.pixel,
  dose.th = 0.02,
  delta.r = 3,
  analysis.th = 0.05,
  local = FALSE,
  local.th = 0.3,
  alias = "",
  description = NULL
)
```

Arguments

vol	"volume" class object, which represents the measured volume.
vol.ref	"volume" class object, which represents the reference volume.
vol.max	Positive number, by default equal to the maximum value of the reference volume. See Details.
dose.th	Positive number, in percent, used to determine the dose difference criterion. See Details.
delta.r	Positive number, in mm. Distance difference criterion.
analysis.th	Positive number, in percent. Only the voxels whose value is greater than or equal analyse.th*vol.max are processed.
local	Boolean. If local = FALSE (default), a global Gamma index is computed, and a local Gamma index otherwise.
local.th	Positive number, in percent. Local threshold, only used if local = TRUE. See Details.
alias	Character string, \$object.alias of the created object.
description	Character string, describing the created object. If description = NULL (default value), it will be set to Gamma index setup.

rt.gamma.index 105

Details

The Gamma index of a voxel n was defined by Low and al [1]. It is computed from the formulae:

$$\gamma_n = min\left(\sqrt{\frac{(D_i - Dref_n)^2}{\Delta D^2} + \frac{{r_i}^2}{\Delta r^2}}\right)$$

whith D_i the measured dose at voxel i, $Dref_n$ the reference dose at voxel n, r_i the distance between voxels i and n, Δr the distance difference criterion equal to delta.r, ΔD the distance difference criterion at voxel n defined as follows:

- If local = FALSE a global Gamma index is computed and $\Delta D = dose.th \cdot vol.max$.
- If local = TRUE, then $\Delta D = dose.th \cdot Dref_n$ when $Dref_n \geq local.th \cdot vol.max$, and $\Delta D = dose.th \cdot local.th \cdot vol.max$ otherwise.

Value

Returns a "volume" class object (see espadon.class for class definitions). The \$vol3D.data field represents the Gamma index. Two fields are added: the \$setup field recalls the calculation setup, and the \$gamma.info field details the rate of Gamma indices below 1, above 1.2 and 1.5, the max and the mean Gamma index.

References

[1] Low DA, Harms WB, Mutic S, Purdy JA (1998). "A technique for the quantitative evaluation of dose distributions." *Medical Physics*, **25**(5), 656–661.

See Also

rt.chi.index

```
# Creation of a reference volume and measured volume
# loading of toy-patient objects (decrease dxyz for better result)
patient <- toy.load.patient (modality = c ("rtdose", "rtstruct"),</pre>
                              roi.name = "ptv", dxyz = c (3, 3, 3))
D.ref <- patient$rtdose[[1]]</pre>
# We will assume that the measured dose is equal to the reference dose shifted
# by one pixel on the x axis
D.meas <- vol.copy (D.ref, alias = "measured_dose")</pre>
D.meas$vol3D.data[1:(D.meas$n.ijk[1] - 1) ,,] <- D.ref$vol3D.data[2:D.ref$n.ijk[1],,]</pre>
gamma <- rt.gamma.index (D.meas, D.ref, delta.r = 6)</pre>
gamma$gamma.info
# Display gamma index at isocenter
G.iso <- patient$rtstruct[[1]]$roi.info$Gz[</pre>
  patient$rtstruct[[1]]$roi.info$name == "ptv"]
display.plane(gamma, view.coord = G.iso,
              bottom.col = c ("#00FF00", "#007F00", "#FF8000", "#FF0000",
                               "#AF0000"),
              bottom.breaks = c (0, 0.8, 1, 1.2, 1.5, gamma$max.pixel),
              bg = "blue", interpolate = FALSE)
```

106 rt.indices.from.bin

 $\begin{tabular}{ll} \it T. indices. from. bin &\it Do simetry, volume, conformity, homogeneity indices from binary selection &\it Conformity for the property of t$

Description

The rt.indices.from.bin function calculates, from a "volume" class object of modality "rtdose", all the standard indicators of radiotherapy, as long as their options are transmitted, for the target and healthy "volume" object of modality "binary".

Usage

```
rt.indices.from.bin(
  vol,
  target.bin.list = NULL,
 healthy.bin.list = NULL,
 T.MAT = NULL,
  presc.dose = NA,
  healthy.tol.dose = NA,
 healthy.weight = 1,
  dosimetry = c("D.min", "D.max", "D.mean", "STD"),
 volume.indices = c("V.tot", "area", "V.prescdose"),
 conformity.indices = c("PITV", "CI.lomax2003", "CN", "NCI", "DSC", "CI.distance",
   "CI.abs_distance", "CDI", "CS3", "ULF", "OHTF", "gCI", "COIN", "G_COSI", "COSI"),
  homogeneity.indices = c("HI.RTOG.max_ref", "HI.RTOG.5_95", "HI.ICRU.max_min",
    "HI.ICRU.2.98_ref", "HI.ICRU.2.98_50", "HI.ICRU.5.95_ref", "HI.mayo2010",
    "HI.heufelder"),
  gradient.indices = c("GI.ratio.50"),
 D.xpc = NULL,
 D.xcc = NULL,
 V.xGy = NULL,
  verbose = TRUE
)
```

Arguments

"volume" class object, of "rtdose" modality. vol target.bin.list list of "volume" class objects, of "binary" modality. The \$object.alias field of each target.bin.list object represents the name of the selected region of interest of the target volume. healthy.bin.list list of "volume" class objects, of "binary" modality. The \$object.alias field of each healthy.bin.list object represents the name of the selected region of interest of the healthy tissues. T.MAT "t.mat" class object, created by load.patient.from.Rdcm or load.T.MAT. If T. MAT = NULL, all \$ref.pseudo of bin.list elements must be equal to vol\$ref.pseudo. presc.dose vector of prescription doses that serve as reference doses for the target RoI. healthy.tol.dose vector of tolerance dose of each healthy RoI.

rt.indices.from.bin 107

```
healthy weight Vector of weight, indicating the importance of the healthy RoI.
dosimetry
                  Vector indicating the requested dose indicators from among 'D.min', 'D.max',
                  'D.mean' and 'STD.' If D. xpc is different from NULL, it will be added.
volume.indices Vector indicating the requested volume indices from among 'V.tot', 'V.prescdose'
                  (i.e. volume over presc.dose) and 'area'. If V.xGy is different from NULL, it
                  will be added.
conformity.indices
                  Vector. Requested conformity indices from among 'PITV', 'PDS', 'CI.lomax2003',
                  'CN', 'NCI', 'DSC', 'CI.distance', 'CI.abs_distance', 'CDI', 'CS3', 'ULF',
                  'OHTF', 'gCI', 'COIN', 'COSI' and 'G_COSI'.
homogeneity.indices
                  Vector. Requested homogeneity indices from among 'HI.RTOG.max_ref', 'HI.RTOG.5_95',
                  'HI.ICRU.max_min', 'HI.ICRU.2.98_ref', 'HI.ICRU.2.98_50', 'HI.ICRU.5.95_ref',
                  'HI.mayo2010' and 'HI.heufelder.'
gradient.indices
                  Vector. Requested gradient indices from among 'GI.ratio.50', 'mGI'.
D.xpc
                  Vector of the percentage of the volume, for which the dose coverage is requested.
                  Vector of the volume in cm^3, for which the dose coverage is requested.
D.xcc
                  Vector of the minimum dose in Gy, received by the volume to be calculated.
V.xGv
                  Boolean. if TRUE (default) a progress bar is displayed.
verbose
```

Value

Return a list of indices dataframe. For details, see rt.indices.from.roi.

See Also

rt.indices.from.roi.

```
# loading of toy-patient objects (decrease dxyz and increase beam.nb for better
# result)
step <- 5
patient <- toy.load.patient (modality = c("rtdose", "rtstruct"), roi.name = "eye",</pre>
                               dxyz = rep (step, 3), beam.nb = 3)
D <- patient$rtdose[[1]]</pre>
struct <- patient$rtstruct[[1]]</pre>
T.MAT <- patient$T.MAT
# creation of the list of target binary volumes
taget.roi.idx <- select.names (struct$roi.info$roi.pseudo, roi.sname = "ptv")</pre>
healthy.roi.idx <- select.names (struct$roi.info$roi.pseudo, roi.sname = "eye")
target.bin.list <- lapply (taget.roi.idx , function (idx) {</pre>
  vr <- nesting.roi (D, struct, roi.idx = idx, xyz.margin = c (5, 5, 5),</pre>
                      T.MAT = T.MAT, alias = struct$roi.info$name[idx])
  b <- bin.from.roi(vr, struct, roi.idx = idx, T.MAT = T.MAT,
                     alias = struct$roi.info$name[idx])
  })
names (target.bin.list) <- struct$roi.info$name[taget.roi.idx]</pre>
healthy.bin.list <- lapply (healthy.roi.idx , function (idx) {</pre>
  vr <- nesting.roi (D, struct, roi.idx = idx, xyz.margin = c (5, 5, 5),</pre>
```

108 rt.indices.from.roi

rt.indices.from.roi Dosimetry, volume, conformity, homogeneity indices from RoI

Description

The rt.indices.from.roi function calculates, from a "volume" class object of modality "rtdose", standard indicators of radiotherapy in relation to the target and healthy RoI, as long as their options are transmitted.

Usage

```
rt.indices.from.roi(
  vol,
  struct = NULL,
 T.MAT = NULL,
  target.roi.name = NULL,
  target.roi.sname = NULL,
  target.roi.idx = NULL,
 healthy.roi.name = NULL,
 healthy.roi.sname = NULL,
 healthy.roi.idx = NULL,
 presc.dose = NA,
 healthy.tol.dose = NA,
 healthy.weight = 1,
 dosimetry = c("D.min", "D.max", "D.mean", "STD"),
  volume.indices = c("V.tot", "area", "V.prescdose"),
  conformity.indices = c("PITV", "PDS", "CI.lomax2003", "CN", "NCI", "DSC",
   "CI.distance", "CI.abs_distance", "CDI", "CS3", "ULF", "OHTF", "gCI", "COIN",
    "G_COSI", "COSI"),
 homogeneity.indices = c("HI.RTOG.max_ref", "HI.RTOG.5_95", "HI.ICRU.max_min",
    "HI.ICRU.2.98_ref", "HI.ICRU.2.98_50", "HI.ICRU.5.95_ref", "HI.mayo2010",
    "HI.heufelder"),
  gradient.indices = c("GI.ratio.50", "mGI"),
 D.xpc = NULL,
 D.xcc = NULL,
 V.xGy = NULL,
  verbose = TRUE
```

Arguments

vol "volume" class object, of "rtdose" modality.

struct "struct" class object.

T.MAT "t.mat" class object, created by load.patient.from.Rdcm or load.T.MAT. If T.MAT

= NULL, struct\$ref.pseudo must be equal to vol\$ref.pseudo.

target.roi.name

Exact name of target RoI in struct object. By default target.roi.name =

NULL. See Details.

target.roi.sname

 $Name\ or\ part\ of\ name\ of\ target\ RoI\ in\ \texttt{struct}\ object.\ By\ default\ \texttt{target.roi.sname}$

= NULL. See Details.

target.roi.idx Value of the index of target RoI that belong to the struct object. By default

target.roi.idx = NULL. See Details.

healthy.roi.name

Exact name of healthy RoI in struct object. By default healthy.roi.name =

NULL.

healthy.roi.sname

 $Name\ or\ part\ of\ name\ of\ healthy\ RoI\ in\ \texttt{struct}\ object.\ By\ default\ \texttt{healthy.roi.sname}$

= NULL.

healthy.roi.idx

Value of the index of healthy RoI that belong to the struct object. By default

healthy.roi.idx = NULL.

presc.dose Vector of prescription doses that serve as reference doses for the target RoI.

healthy.tol.dose

Vector of tolerance doses of each healthy RoI.

healthy.weight Vector of weights, indicating the importance of the healthy RoI.

dosimetry Vector indicating the requested dose indicators from among 'D.min', 'D.max',

'D.mean' and 'STD.' If D.xpc is different from NULL, it will be added.

volume.indices Vector indicating the requested volume indices from among 'V.tot', 'V.prescdose'

(i.e. volume over presc.dose) and 'area'. If V.xGy is different from NULL, it

will be added.

conformity.indices

Vector. Requested conformity indices from among 'PITV', 'PDS', 'CI.lomax2003',

'CN', 'NCI', 'DSC', 'CI.distance', 'CI.abs_distance', 'CDI', 'CS3', 'ULF',

'OHTF', 'gCI', 'COIN', 'COSI' and 'G_COSI'.

homogeneity.indices

Vector. Requested homogeneity indices from among 'HI.RTOG.max_ref', 'HI.RTOG.5_95',

'HI.ICRU.max_min', 'HI.ICRU.2.98_ref', 'HI.ICRU.2.98_50', 'HI.ICRU.5.95_ref',

'HI.mayo2010' and 'HI.heufelder.'

gradient.indices

Vector. Requested gradient indices from among 'GI.ratio.50', 'mGI'.

D. xpc Vector of the percentage of the volume, for which the dose coverage is requested.

D. xcc Vector of the volume in cm^3 , for which the dose coverage is requested.

V.xGy Vector of the minimum dose in Gy, received by the volume to be calculated.

verbose Boolean. if TRUE (default) a progress bar is displayed.

Details

If target.roi.name, target.roi.sname, and target.roi.idx are all set to NULL, all RoI containing 'ptv' (if they exist) are selected.

If target.roi.name, target.roi.sname, and target.roi.idx are all set to NULL, all RoI that include 'ptv', if they exist, are selected by default.

If healthy.roi.name, healthy.roi.sname, and healthy.roi.idx are all set to NULL, no healthy RoI are selected.

Value

Return a list containing (if requested)

- -dosimetry: dataframe containing, for all target and healthy structures:
 - the requested dosimetry: D.min (Gy), D.max (Gy), D.mean (Gy) and STD (Gy), respectively the minimum, maximum, mean and standard deviation of the dose in the regions of interest.
 - the requested \$D.x%: (Gy) Dose covering x percent of structure volume.
 - the requested D.xcc: (Gy) Dose covering x (cm^3) of structure volume.
- -volume: dataframe containing, for all target and healthy structures, and isodoses:
 - the requested volume.indices: V_tot (cm^3) (except for isodose) the total volume of the regions of interest, area (cm^2) (except for isodose) their surface areas, V.prescdose (cm^3) the volumes receiving at least presc.dose Gy,
 - the requested $V.xGy(cm^3)$: volumes receiving at least x.Gy.
- -conformity: dataframe containing, if requested,
 - PITV: (1) Prescription Isodose Target Volume, or conformity index defined by E.Shaw [1]

$$PITV = \frac{V_{presc.dose}}{V_{target}}$$

• PDS: (1) Prescription Dose Spillage defined by SABR UK Consortium 2019 [2]

$$PDS = \frac{V_{presc.dose}}{V_{target \geq presc.dose}} = \frac{V_{presc.dose}}{V_{target} \ \cap \ V_{presc.dose}}$$

• CI.lomax2003: (1) Conformity Index defined by Lomax and al [3]

$$CI_{lomax2003} = \frac{V_{target \ \geq \ presc.dose}}{V_{presc.dose}} = \frac{V_{target \ \cap \ }V_{presc.dose}}{V_{presc.dose}}$$

• CN: (1) Conformation Number defined by *Van't Riet and al* [4]. It corresponds to conformity index defined by *Paddick* [5]

$$CN = CI_{paddick2000} = \frac{V_{target \geq presc.dose}^2}{V_{target} \cdot V_{presc.dose}} = \frac{(V_{target} \cap V_{presc.dose})^2}{V_{target} \cdot V_{presc.dose}}$$

• NCI: (1) New conformity index, inverse of CN, defined by Nakamura and al [6]

$$NCI = \frac{1}{CN}$$

• DSC: (1) Dice Similarity Coefficient [7]

$$DSC = 2 \cdot \frac{V_{target \geq presc.dose}}{V_{target} + V_{presc.dose}} = 2 \cdot \frac{V_{target} \cap V_{presc.dose}}{V_{target} + V_{presc.dose}}$$

• CI. distance: (1) Conformity Index based on distance defined by Park and al [8]

$$CI.distance = \frac{100}{N} \sum_{i=1}^{N} \frac{dist_{S_{presc.dose} \to G_{target}} - dist_{S_{target} \to G_{target}}}{dist_{S_{target} \to G_{target}}}$$

where $dist_{S_{presc.dose} \to G_{target}}$ is the distance between the surface of the prescription dose volume and the centroid of the target, and $dist_{S_{target} \to G_{target}}$ the surface of the target volume and the centroid of the target. N is the number of directions where the distances are calculated. These directions are computed every 1° . If the centroid is not within the target volume, then CI.distance = NA.

• CI.abs_distance: (1) Conformity Index based on distance defined by Park and al [8]

$$CI.abs_distance = \frac{100}{N} \sum_{i=1}^{N} \frac{|dist_{S_{presc.dose}} \to G_{target} - dist_{S_{target}} \to G_{target}|}{dist_{S_{target}} \to G_{target}}$$

• CDI: (1) Conformity Distance Index defined by Wu and al [9]

$$CDI = 2 \frac{V_{presc.dose} + V_{target} - 2 \ V_{target} \ge presc.dose}{S_{target} + S_{presc.dose}} = \frac{V_{presc.dose} + V_{target} - 2 \ \cdot \ V_{target} \ \cap \ V_{presc.dose}}{S_{target} + S_{presc.dose}}$$

where S_{target} is the surface of the target volume and $S_{presc.dose}$ is the surface of the prescription dose volume.

• CS3: (1) Triple Point Conformity Scale defined by Ansari and al [10]

$$CS3 = \frac{V_{0.95 \cdot presc.dose} + V_{presc.dose} + V_{1.05 \cdot presc.dose}}{3 \cdot V_{target}}$$

• ULF: (1) Underdosed lesion factor defined by Lefkopoulos and al [11]

$$ULF = \frac{V_{target} < presc.dose}{V_{target}} = \frac{V_{target} \cap \overline{V_{presc.dose}}}{V_{target}}$$

• OHTF:(1) Overdosed healthy tissues factor defined by Lefkopoulos and al [11]

$$OHTF = \frac{\sum V_{healthy \ge presc.dose}}{V_{target}} = \frac{\sum V_{healthy} \ \cap \ V_{presc.dose}}{V_{target}}$$

• gCI: (1) Geometric Conformity Index defined by Lefkopoulos and al [11]

$$qCI = ULF + OHTF$$

• COIN: Conformity Index defined by Baltas and al [12]

$$COIN = \frac{V_{target \geq presc.dose}^2}{V_{target} \cdot V_{presc.dose}} \cdot \prod^{N_{healthy}} \left(1 - \frac{V_{healthy \geq presc.dose}}{V_{healthy}}\right)$$

• gCOSI: generalized COSI, defined by Menhel and al [13].

$$gCOSI = 1 - \sum^{N_{healthy}} healthy.weight \cdot \frac{\frac{V_{healthy} \geq healthy.tol.dose}{V_{healthy}}}{\frac{V_{target} \geq presc.dose}{V_{target}}}$$

- COSI: if "COSI" is requested in conformity.indices, it returns a dataframe of Critical Organ Scoring Index for each healthy organ, at each presc.dose, and for each target. COSI is defined by Menhel and al [13]

$$COSI = 1 - \frac{\frac{V_{healthy} \ge healthy.tol.dose}{V_{healthy}}}{\frac{V_{target} \ge presc.dose}{V_{target}}}$$

- homogeneity: dataframe containing
 - HI.RTOG.max_ref: (1) Homogeneity Index from RTOG defined by E.Shaw [1]

$$HI.RTOG.max_ref = \frac{D_{max}}{presc.dose}$$

where D_{max} is the maximum dose in the target volume.

• HI.RTOG.5_95: (1) Homogeneity Index from RTOG [1]

$$HI.RTOG.5_95 = \frac{D.5pc}{D.95pc}$$

where D.5pc and D.95pc are respectively the doses at 5% and 95% of the target volume in cumulative dose-volume histogram.

• HI.ICRU.max_min: (1) Homogeneity Index from ICRU report 62 [14]

$$HI.ICRU.max_min = \frac{D_{max}}{D_{min}}$$

where D_{max} and D_{min} are respectively the maximum and the minimum dose in the target volume.

• HI.ICRU.2.98_ref: (1) Homogeneity Index from ICRU report 83 [15]

$$HI.ICRU.2.98_ref = 100 \frac{D.2pc - D.98pc}{presc.dose}$$

where D.2pc and D.98pc are respectively the doses at 2% and 98% of the target volume in cumulative dose-volume histogram.

• HI.ICRU.2.98_50: (1) Homogeneity Index from ICRU report 83 [15]

$$HI.ICRU.2.98_50 = 100 \frac{D.2pc - D.98pc}{D.50pc}$$

where D.2pc, D.98pc and D.50pc are respectively the doses at 2%, 98% and 50% of the target volume in cumulative dose-volume histogram.

• HI.ICRU.5.95_ref: (1) Homogeneity Index from ICRU report 83 [15]

$$HI.ICRU.5.95_ref = 100 \frac{D.5pc - D.95pc}{presc.dose}$$

where D.5pc and D.95pc are respectively the doses at 5% and 95% of the target volume in cumulative dose-volume histogram.

• HI. mayo2010: (1) Homogeneity Index defined by Mayo and al [16]

$$HI.mayo2010 = \sqrt{\frac{D_{max}}{presc.dose} \, \cdot \, (1 + \frac{\sigma_D}{presc.dose})}$$

where D_{max} is the maximum dose in the target volume, and σ_D the standard deviation of the dose in the target volume.

• HI. heufelder: (1) Homogeneity Index defined by *Heufelder and al* [17]

$$HI.heufelder = e^{-0.01 \cdot (1 - \frac{\mu_D}{presc.dose})^2} \cdot e^{-0.01 \cdot (\frac{\sigma_D}{presc.dose})^2}$$

where μ_D and σ_D are respectively the mean and the standard deviation of the dose in the target volume.

- gradient : dataframe containing
 - GI.ratio.50: Gradient Index based on volumes ratio defined by Paddick and Lippitz [18]

$$GI.ratio.50 = \frac{V_{0.5 \cdot presc.dose}}{V_{presc.dose}}$$

• mGI: Modified Gradient Index defined by SABR UK Consortium 2019 [2]

$$mGI = \frac{V_{0.5 \ \cdot \ presc.dose}}{V_{target \ \geq \ presc.dose}} = \frac{V_{0.5 \ \cdot \ presc.dose}}{V_{target \ \cap \ V_{presc.dose}}}$$

References

- [1] Shaw E, Kline R, Gillin M, Souhami L, Hirschfeld A, Dinapoli R, Martin L (1993). "Radiation therapy oncology group: Radiosurgery quality assurance guidelines." *International journal of radiation oncology, biology, physics*, **27**(5), 1231-1239. ISSN 0360-3016, doi:10.1016/0360-3016(93)90548A.
- [2] UK SABR Consortium (Online; accessed 2022-04-01). "Stereotactic Ablative Radiation Therapy (SABR): a resource. v6.1, January 2019." https://www.sabr.org.uk/wp-content/uploads/2019/04/SABRconsortium-guidelines-2019-v6.1.0.pdf.
- [3] Lomax NJ, Scheib SG (2003). "Quantifying the degree of conformity in radiosurgery treatment planning." *International journal of radiation oncology, biology, physics*, **55**(5), 1409-1419. ISSN 0360-3016, doi:10.1016/S03603016(02)045996.
- [4] Riet AV, Mak AC, Moerland MA, Elders LH, Van der Zee W (1997). "A conformation number to quantify the degree of conformality in brachytherapy and external beam irradiation: Application to the prostate." *International journal of radiation oncology, biology, physics*, **37**(3), 731-736. ISSN 0360-3016, doi:10.1016/S03603016(96)006013.
- [5] Paddick I (2000). "A simple scoring ratio to index the conformity of radiosurgical treatment plans. Technical note." *Journal of neurosurgery*, **93 Suppl 3**, 219-222.
- [6] Nakamura J, Verhey L, Smith V, Petti P, Lamborn K, Larson D, Wara W, Mcdermott M, Sneed P (2002). "Dose conformity of Gamma Knife radiosurgery and risk factors for complications." *International journal of radiation oncology, biology, physics*, **51**, 1313-9. doi:10.1016/S0360-3016(01)017576.
- [7] Dice LR (1945). "Measures of the Amount of Ecologic Association Between Species." *Ecology*, **26**(3), 297–302. ISSN 00129658, 19399170.
- [8] Park JM, Park S, Ye S, Kim J, Carlson J, Wu H (2014). "New conformity indices based on the calculation of distances between the target volume and the volume of reference isodose." *The British journal of radiology*, **87**, 20140342. doi:10.1259/bjr.20140342.
- [9] Wu Q, Wessels BW, Einstein DB, Maciunas RJ, Kim EY, Kinsella TJ (2003). "Quality of coverage: Conformity measures for stereotactic radiosurgery." *Journal of Applied Clinical Medical Physics*, **4**, 374-381.
- [10] Ansari S, Satpathy S, Singh P, Lad S, Thappa N, Singh B (2018). "A new index: Triple Point Conformity Scale (CS3) and its implication in qualitative evaluation of radiotherapy plan." *Journal of Radiotherapy in Practice*, 17, 1-4. doi:10.1017/S1460396917000772.

[11] Lefkopoulos D, Dejean C, balaa ZE, Platoni K, Grandjean P, Foulquier J, Schlienger M (2000). "Determination of dose-volumes parameters to characterise the conformity of stereotactic treatment plans." In chapter XIII, 356-358. Springer Berlin Heidelberg. ISBN 978-3-540-67176-3, doi:10.1007/9783642597589_135.

- [12] Baltas D, Kolotas C, Geramani KN, Mould RF, Ioannidis G, Kekchidi M, Zamboglou N (1998). "A conformal index (COIN) to evaluate implant quality and dose specification in brachytherapy." *International journal of radiation oncology, biology, physics*, **40 2**, 515-24. doi:10.1016/s0360-3016(97)007323.
- [13] Menhel J, Levin D, Alezra D, Symon Z, Pfeffer R (2006). "Assessing the quality of conformal treatment planning: a new tool for quantitative comparison." *Physics in Medicine and Biology*, 51(20), 5363–5375.
- [14] Landberg T, Chavaudra J, Dobbs J, Gerard J, Hanks G, Horiot J, Johansson K, Möller T, Purdy J, Suntharalingam N, Svensson H (1999). "ICRU Report 62: Prescribing, Recording and Reporting Photon Beam Therapy (Supplement to ICRU Report 50),3. Absorbed Doses." *Reports of the International Commission on Radiation Units and Measurements*, os-32(1), 21-25.
- [15] ICRU (2010). "Report 83: Prescribing, Recording, and Reporting Photon-Beam Intensity-Modulated Radiation Therapy (IMRT)." *Reports of the International Commission on Radiation Units and Measurements*, **10**(1), 1-3.
- [16] Mayo CS, Ding L, Addesa A, Kadish S, Fitzgerald TJ, Moser R (2010). "Initial Experience With Volumetric IMRT (RapidArc) for Intracranial Stereotactic Radiosurgery." *International Journal of Radiation Oncology*Biology*Physics*, **78**(5), 1457-1466. ISSN 0360-3016, doi:10.1016/j.ijrobp.2009.10.005.
- [17] Heufelder J, Zink K, Scholz M, Kramer K, Welker K (2003). "Eine Methode zur automatisierten Bewertung von CT-basierten Bestrahlungsplanen in der perkutanen Strahlentherapie." *Zeitschrift fur Medizinische Physik*, **13**(4), 231-239. ISSN 0939-3889, doi:10.1078/0939388900175.
- [18] Paddick I, Lippitz BE (2006). "A simple dose gradient measurement tool to complement the conformity index." *Journal of neurosurgery*, **105 Suppl**, 194-201.

All this references are compiled by

- Kaplan LP, Korreman SS (2021). "A systematically compiled set of quantitative metrics to describe spatial characteristics of radiotherapy dose distributions and aid in treatment planning." *Physica Medica*, **90**, 164-175. ISSN 1120-1797, doi:10.1016/j.ejmp.2021.09.014. and
- Patel G, Mandal A, Choudhary S, Mishra R, Shende R (2020). "Plan evaluation indices: A journey of evolution." *Reports of Practical Oncology & Radiotherapy*, **25**. doi:10.1016/j.rpor.2020.03.002..

See Also

rt.indices.from.bin.

save.T.MAT

verbose = FALSE)

indices

save.T.MAT

Save a T.MAT class object

Description

The save.T.MAT function saves the data required by load.T.MAT, load.patient.from.dicom or load.patient.from.Rdcm to generate T.MAT, as pre-formatted Rdcm files.

Usage

```
save.T.MAT(T.MAT, dirname)
```

Arguments

T. MAT "t.mat" class object to save.

dirname Directory where new reg .Rdcm files will be saved.

Details

Reg files from DICOM files cannot be updated with the save. T. MAT function. Only transfer matrices added with ref.add or ref.cutplane.add will be saved.

Value

Returns TRUE, if all reg files generating T.MAT are saved.

116 save.to.Rdcm

save.to.Rdcm

Save a espadon object in a pre-formatted *.Rdcm file

Description

The function save.to.Rdcm allows you to save an object created by **espadon** in a pre-formatted *.Rdcm file. This object will also be accessible by the load.patient.from.Rdcm function.

Usage

```
save.to.Rdcm(obj, object.name = obj$object.alias, dirname = obj$file.dirname)
```

Arguments

```
obj espadon object of class "volume", "struct", "mesh", "histo", "dvh", "histo2D".

object.name Character string, representing the name of the object, default to obj$object.alias.

dirname Directory where new files from obj will be saved.
```

Value

Returns TRUE, if paste0(object.name, ".Rdcm") exists in dirname.

Returns FALSE, if object.name is not a valid file name, or if the file that is created would replace a *.Rdcm file created by dicom.to.Rdcm.converter.

Note

save.to.Rdcm can not replace an *.Rdcm file created by dicom.to.Rdcm.converter.

select.names 117

select.names	Regions of Interest (RoI) indices

Description

The select.names function allows you to select words from a vector of words, according to several criteria, eliminating spaces and case.

Usage

```
select.names(names, roi.name = NULL, roi.sname = NULL, roi.idx = NULL)
```

Arguments

names	Words vector
roi.name	Vector of words to compare to names. By default roi.name = NULL. See Details
roi.sname	Vector of words or parts of words to compare. By default roi.sname = NULL. See Details
roi.idx	Index vector. By default roi.idx = NULL. See Details.

Details

If roi.name, roi.sname, and roi.idx are all NULL, then all RoI are selected.

Value

Returns the indices of the elements of the word vector names satisfying one or more of the following conditions:

- ASCII // TRANSLIT transcriptions, without spaces, of names and roi.name, are identical.
- ASCII // TRANSLIT transcriptions, without spaces of roi.sname are identical to part of ASCII // TRANSLIT transcriptions, without spaces of names.
- names indices belong to the index vector roi.idx.

```
# loading patient objects
names <- c ("Eye left", "EyeR", "OPTICAL nerve L", "opical nervR", "chiasma")
# RoI selection.
select.names (names = names, roi.name = c("eye left", "eye right"))
select.names (names = names, roi.sname = c("eye", "ner"))
select.names (names = names, roi.idx = 4:9)</pre>
```

118 struct.clustering

struct.clustering Clustering volumes by RoI

Description

The struct.clustering function creates a new volume in which voxels are clustered and labeled by region of interest defined in an rt-struct.

Usage

```
struct.clustering(
  vol,
  struct,
  roi.name = NULL,
  roi.sname = NULL,
  roi.idx = NULL,
  T.MAT = NULL,
  alias = "",
  description = NULL,
  verbose = TRUE
)
```

Arguments

vol	"volume" class object.
struct	"struct" class object.
roi.name	Vector of exact names of the RoI in the struct object. By default roi.name = NULL. See Details.
roi.sname	Vector of names or parts of names of the RoI in the struct object. By default roi.sname = NULL. See Details.
roi.idx	Vector of indices of the RoI that belong to the struct object. By default roi.idx = NULL. See Details.
T.MAT	"t.mat" class object, created by load.patient.from.Rdcm or load.T.MAT. If T.MAT = NULL, struct\$ref.pseudo must be equal to vol\$ref.pseudo.
alias	Character string, \$alias of the created object.
description	Character string, describing the created object. If description = NULL (default value), it will be set to paste (struct\$object.alias, "clustering")
verbose	Boolean. if TRUE (default), the RoI studied are listed.

Details

If roi.name, roi.sname, and roi.idx are all set to NULL, all RoI are selected.

Value

Returns a "volume" class object (see espadon.class for class definitions), of "cluster" modality. This object contains the \$cluster.info field, detailing the label and volumes in cm^3 of the different clusters. Note that the label NA or value 0 is used for the voxels which are not contained in any RoI (air for instance).

struct.from.bin

See Also

get.roi.connection

Examples

struct.from.bin

Creation of struct class object from a binary volume

Description

The struct.from.bin function creates a struct object with a unique RoI, defined by the contours of binary volume.

Usage

```
struct.from.bin(
  vol,
  roi.name = vol$description,
  roi.nb = 1,
  roi.color = "#379DA2",
  roi.type = c("", "EXTERNAL", "PTV", "CTV", "GTV", "TREATED_VOLUME", "IRRAD_VOLUME",
        "OAR", "BOLUS", "AVOIDANCE", "ORGAN", "MARKER", "REGISTRATION", "ISOCENTER",
        "CONTRAST_AGENT", "CAVITY", "BRACHY_CHANNEL", "BRACHY_ACCESSORY", "BRACHY_SRC_APP",
        "BRACHY_CHNL_SHLD", "SUPPORT", "FIXATION", "DOSE_REGION", "CONTROL",
        "DOSE_MEASUREMENT"),
    external.only = FALSE,
    alias = "",
    description = paste("RoI from", vol$object.alias)
)
```

```
vol "volume" class object, of binary modality.
roi.name Character string, representing the name of created RoI.
```

120 struct.in.new.ref

roi.nb	Positive integer, representing the number of created RoI.
roi.color	Color of the created RoI, in hex code format ("#RRGGBB").
roi.type	Type of RoI, from among "", "EXTERNAL", "PTV", "CTV", "GTV", "TREATED_VOLUME", "IRRAD_VOLUME", "OAR", "BOLUS", "AVOIDANCE", "ORGAN", "MARKER", "REGISTRATION", "ISOCENTER", "CONTRAST_AGENT", "CAVITY", "BRACHY_CHANNEL "BRACHY_ACCESSORY", "BRACHY_SRC_APP", "BRACHY_CHNL_SHLD", "SUPPORT", "FIXATION", "DOSE_REGION", "CONTROL" and "DOSE_MEASUREMENT"
external.only	Boolean. If TRUE, only external contours are kept.
alias	Character string, \$alias of the created object.
description	Character string, describing the created object.

Value

Returns a "struct" class object (see espadon.class for class definition), including the unique roi. name as region of interest.

Examples

struct.in.new.ref Change of frame of reference of a "struct" class object.

Description

The struct.in.new.ref function allows you to change the frame of reference of a struct.

Usage

```
struct.in.new.ref(struct, new.ref.pseudo, T.MAT, alias = "")
```

struct	"struct" class object.
new.ref.pseudo	pseudonym of the frame of reference in which the struct should be located. This new.ref.pseudo must exist in the T.MAT list.
T.MAT	"t.mat" class object, created by load.patient.from.dicom, load.patient.from.Rdcm, load.T.MAT or ref.add.
alias	Character string, \$alias of the created object.

struct.merge 121

Value

Returns "struct" class object in the new frame of reference new.ref.pseudo.

See Also

vol.in.new.ref

Examples

struct.merge

Merging of structures into a new structure

Description

The struct.merge function merges two structures into a new one. It is useful for comparing contours, for example.

Usage

```
struct.merge(
  ref.struct,
  add.struct,
  roi.name = NULL,
  roi.sname = NULL,
  roi.idx = NULL,
  suffix = "",
  alias = "",
  description = ""
)
```

ref.struct	struct class object. All RoI of this structure are kept.
add.struct	struct class object. Only the selected RoI are kept for merging.
roi.name	Vector of exact names of the RoI in the add.struct object. By default roi.name = NULL. See Details.
roi.sname	Vector of names or parts of names of the RoI in the add.struct object. By default roi.sname = NULL. See Details.
roi.idx	Vector of indices of the RoI that belong to the add.struct object. By default roi.idx = NULL. See Details.
suffix	Character string. '-suffix' is added to RoI name.
alias	Character string, \$alias of the resulted object.
description	Character string, describing the the resulted object.

122 toy.dicom.raw

Details

If roi.name, roi.sname, and roi.idx are all NULL, then all RoI of add.struct are selected.

Value

Returns a struct class object. See espadon.class for class definitions.

Note

Beware that, when merging structures, some RoI may have same name or roi.info\$roi.pseudo. In this case struct.merge prints a warning message. Consider changing suffix to avoid the ambiguity.

See Also

struct.from.bin.

Examples

```
# loading of toy-patient objects (decrease dxyz and increase beam.nb for
# better result)
step <- 5
patient <- toy.load.patient (modality = c("rtdose"),</pre>
                              dxyz = rep (step, 3), beam.nb = 3)
D <- patient$rtdose[[1]]</pre>
# isodose 50% Dmax Gy and 90% Dmax
bin50 <- bin.from.vol (D, min = 0.5 * D$max.pixel)</pre>
bin90 <- bin.from.vol (D, min = 0.9 * D$max.pixel)</pre>
S.isodose50 <- struct.from.bin (bin50, roi.name = "50pc" ,</pre>
                                 roi.color = "#00FFFF")
S.isodose90 <- struct.from.bin (bin90, roi.name = "90pc"</pre>
                                  roi.color = "#FFFF00")
S <- struct.merge (S.isodose50, S.isodose90, alias = "isodose",
                    description = paste ("isodose of", D$object.alias))
# Dmax location :
z.dmax <- get.xyz.from.index(which (D$vol3D.data == D$max.pixel), D)[1,3]</pre>
display.plane(top = D, struct = S, view.coord = z.dmax, legend.shift = -50)
```

toy.dicom.raw

toy DICOM raw data

Description

The toy.dicom.raw loads raw data from a dummy DICOM file. It is used for the test.

Usage

```
toy.dicom.raw()
```

Value

Returns the raw data of a dummy DICOM file of rtplan modality.

toy.load.patient 123

Examples

```
toy.dicom.raw ()
```

toy.load.patient

Load a toy patient for test

Description

The toy.load.patient creates a dummy "patient" class object. It is used for the test.

Usage

```
toy.load.patient(
  modality = c("ct", "mr", "rtdose", "rtstruct"),
  roi.name = c("eye", "optical nerve", "brain", "labyrinth processing unit",
        "energy unit", "gizzard", "ghost container", "exhaust valve"),
  dxyz = c(1, 1, 1),
  beam.nb = 7
)
```

Arguments

modality	String vector, whose elements are chosen among the modalities "ct", "mr", "rt-struct" and "rtdose".
roi.name	String vector, whose elements are chosen among the regions of interest (RoI) "eye", "optical nerve", "brain", "labyrinth processing unit", "energy unit", "gizzard", "ghost container" and "exhaust valve". Note that the RoI "couch", "patient" and "ptv" are still present.
dxyz	Vector of length 3, representing the x, y, z steps in mm, between ct, mr and rtdose voxels.
beam.nb	Positive integer. Number of radiotherapy beams in rtdose modality.

Value

Returns an toy object of "patient" class, containing the modalities defined in modality. See espadon.class for class definitions.

```
# loading of toy-patient objects (decrease dxyz for better result)
step <- 5
pat <- toy.load.patient (dxyz = rep (step, 3), beam.nb = 2)
str (pat, max.level = 2)</pre>
```

124 vol.copy

vector.product Vector product of two vector	vector.product	Vector product of two vectors
---	----------------	-------------------------------

Description

Vector product of two vectors

Usage

```
vector.product(v1, v2)
```

Arguments

v1	Vector of x, y, z coordinates
v2	Vector of x, y, z coordinates

Value

Returns the x, y, z coordinates of the vector product of v1 and v2

have the number of vol.

Examples

```
vector.product(c\ (1,\ \emptyset,\ \emptyset),\ c\ (\emptyset,\ 1,\ \emptyset))
```

vol.copy

Creating a volume from another one

Description

The vol.copy function creates a "volume" class object, with the same grid as the vol volume object.

Usage

```
vol.copy(vol, alias = "", modality = NULL, description = NULL, number = NULL)
```

vol	"volume" class object, template of the created object.
alias	Character string, \$object.alias of the created object.
modality	Character string, modality of the created volume. If modality = NULL, then the created object will have the modality of vol.
description	Character string, description of the returned object. If descritption = NULL, then the created object will have the description of vol.
number	number of the returned volume. If number = NULL, then the returned object will

vol.create 125

Value

Returns a "volume" class object (see espadon.class for class definitions), with the same grid as vol, in which \$vol3D.data is initialized to NA.

Examples

```
# loading of toy-patient objects (decrease dxyz for better result)
patient <- toy.load.patient (modality = "ct", roi.name = "",dxyz = c (4, 4, 4))
CT <- patient$ct[[1]]

# creating a volume
vol.from.CT <- vol.copy (CT, alias = "ct reference")
str (vol.from.CT)</pre>
```

vol.create

Volume creating

Description

The vol.create function creates a volume object from a user-defined grid.

Usage

```
vol.create(
    n.ijk,
    dxyz,
    mid.pt = NULL,
    pt000 = NULL,
    default.value = NA,
    ref.pseudo = "",
    frame.of.reference = "",
    alias = "",
    modality = "",
    description = "",
    number = 0
)
```

Arguments

n.ijk	Vector of length 3, representing the number of elements on the i, j and k axes.	
dxyz	Vector of length 3, representing the x, y, z steps in mm, between voxels. See details.	
mid.pt	Vector of length 3, representing the x, y, z coordinates of the midpoint of the volume. See details.	
pt000	Vector of length 3, representing the x, y, z coordinates of the first voxel of the first plane.	
default.value	Numerical or boolean value, representing the default value of the voxels.	
ref.pseudo	Character string, frame of reference pseudonym of the created object.	
frame.of.reference		

Character string, frame of reference of the created object.

126 vol.from.bin

alias	Character string, \$alias of the created object.
modality	Character string, \$modality of the created object.
description	Character string, describing the the created object.
number	Integer, by default set to 0, number of the created object.

Details

If mid.pt and pt000 are both equal to NULL, then mid.pt = c(0, 0, 0) by default. If mid.pt and pt000 are both different from NULL, then only mid.pt is taken into account.

Value

Returns a "volume" class object (see espadon.class for class definitions), in which the grid is defined by pt000 or mid.pt, dxyz n.ijk. If default.value are initialized to FALSE, then modality = "binary". The orientation of the patient is orthonormal to the grid.

Examples

vol.from.bin

Volume class object according to binary selection

Description

The vol.from.bin function selects a part of a "volume" class object of "binary" modality which has the same grid. It is especially useful to restrict voxel data in region of interest.

Usage

```
vol.from.bin(vol, sel.bin, alias = "", description = NULL)
```

Arguments

vol "volume" class object, containing data to restrict.

sel.bin "volume" class object, of "binary" modality. vol and sel.bin must have the

same grid.

alias Character string, \$alias of the created object

description Character string, describing the created object. If description = NULL (default

value), it will be set to paste (vol\$object.alias, "from", sel.bin\$object.alias)

Value

Returns a "volume" class object (see espadon.class for class definitions), in which non-selected voxels have the value NA, and selected voxels have the original value of vol.

vol.gradient 127

Examples

vol.gradient

Gradient of a volume

Description

The vol.gradient function calculates the 3D gradient of a "volume" class object

Usage

```
vol.gradient(vol, alias = "", description = NULL)
```

Arguments

vol "volume" class object.

alias Character string, \$alias of the created object.

description Character string, describing the created object. If description = NULL (default

value), it will be set to paste (vol\$object.alias, "gradient").

Value

Returns a "volume" class object (see espadon.class for class definitions), with the same grid and modality as vol, gradient of vol.

128 vol.median

vol.in.new.ref

Change of frame of reference of a volume

Description

The vol.in.new.ref function allows you to change the frame of reference of a volume.

Usage

```
vol.in.new.ref(vol, new.ref.pseudo, T.MAT, alias = "", description = NULL)
```

Arguments

vol "volume" class object.

new.ref.pseudo pseudonym of the frame of reference in which the volume should be located. This new.ref.pseudo must exist in the T.MAT list.

T.MAT "t.mat" class object, created by load.patient.from.dicom, load.patient.from.Rdcm, load.T.MAT or ref.add.

alias Character string, \$alias of the created object.

description Character string, describing the created object. If description = NULL (default

value), it will be that of the vol volume.

Value

Returns "volume" class object in the new frame of reference new.ref.pseudo.

See Also

struct.in.new.ref

Examples

vol.median

Median filter on a volume

Description

The vol.median function applies a 26-connectivity median filter on all the voxels of a "volume" class object.

Usage

```
vol.median(vol, alias = "", description = NULL)
```

vol.regrid 129

Arguments

vol "volume" class object.

alias Character string, \$alias of the created object.

description Character string, describing the created object. If description = NULL (default

value), it will be set to paste (vol\$object.alias, "median").

Value

Returns a "volume" class object (see espadon.class for class definitions), with the same grid and modality as vol, in which voxels are filtered by a 26-connectivity median filter.

Examples

vol.regrid

Transform the grid of a volume class object into the grid of another

Description

The vol. regrid function transforms the grid of a volume according to the grid of another.

Usage

```
vol.regrid(
  vol,
  back.vol,
  T.MAT = NULL,
  interpolate = TRUE,
  alias = "",
  description = NULL,
  verbose = TRUE
```

Arguments

vol "volume" class object to regrid.

back.vol "volume" class object whose grid will be used for regriding. Its \$ref.pseudo

must exist in the T. MAT list.

T.MAT "t.mat" class object, created by load.patient.from.dicom, load.patient.from.Rdcm,

load.T.MAT or ref.add. If T.MAT = NULL, back.vol\$ref.pseudo must be equal

to vol\$ref.pseudo.

130 vol.sum

interpolate Boolean, default to TRUE. If interpolate = TRUE, a trilinear interpolation of the value of the voxels, relative to the values of adjacent voxels, is performed.

Character string, \$alias of the created object.

Character string, describing the treated object. If description = NULL, it will be that of vol.

Verbose Boolean. if TRUE (default) a progress bar is displayed.

Value

Returns a copy of vol, in which grid is that of back.vol.

Examples

```
# loading of toy-patient objects (decrease dxyz and increase beam.nb for
# better result)
step <- 5
patient <- toy.load.patient (modality = c ("mr", "rtdose"),</pre>
                              dxyz = rep (step, 3), beam.nb = 4)
MR <- patient$mr[[1]]</pre>
D <- patient$rtdose[[1]]</pre>
# change grid
D.on.MR <- vol.regrid (vol = D, back.vol = MR, interpolate = TRUE,
                       T.MAT = patient$T.MAT, alias = "",
                        description = NULL, verbose = FALSE)
# maximum dose location
max.dose.in.MR <- get.xyz.from.index (which.max (D.on.MR$vol3D.data), D.on.MR)</pre>
display.plane (bottom = MR, view.coord = max.dose.in.MR[3],
               top= D.on.MR, bottom.col = grey.colors(255, start = 0, end = 1),
               bottom.breaks = seq (0, 500, length.out = 256),
               bg = "#00ff00", interpolate = FALSE)
```

vol.sum

Sum of 2 volumes

Description

The vol. sum function adds two "volume" class objects of the same grid and of the same modality.

Usage

```
vol.sum(vol1, vol2, alias = "", description = NULL)
```

vol1, vol2	"volume" class objects. The 2 volumes must have the same modality, and the same grid (i.e. share the same position of the voxels).
alias	Character string, \$alias of the created object.
description	Character string, describing the created object. If description = NULL (default value), it will be set to paste (vol1\$object.alias, "+", vol2\$object.alias).

xlsx.from.dcm 131

Value

Returns a "volume" class object (see espadon.class for class definitions), with the same grid and modality as vol1 and vol2, sum of vol1 and vol2.

Examples

xlsx.from.dcm

Converting DICOM files to .xlsx files

Description

The xlsx.from.dcm function creates an Excel file from DICOM files.

Usage

```
xlsx.from.dcm(
  dcm.filenames,
  xlsx.filenames,
  multipage = TRUE,
  txt.sep = "\\",
  txt.length = 100,
  tag.dictionary = dicom.tag.dictionary()
)
```

Arguments

dcm.filenames	String vector, representing the list of full names of DICOM files.
xlsx.filenames	String vector, representing the list of full names of created *.xlsx files. If multipage = TRUE, only the xlsx.filenames[1] is used.
multipage	Boolean. If TRUE, all dcm. filenames are converted into multiple pages of the same *.xlsx file.
txt.sep	String. Used if as.txt = TRUE. Separator of the tag value elements.
txt.length	Positive integer. Used if as.txt = TRUE. Maximum number of letters in the representation of the TAG value.
tag.dictionary	Dataframe, by default equal to dicom.tag.dictionary, whose structure it must keep. This dataframe is used to parse DICOM files.

Value

Returns a boolean vector, establishing the existence of the created Excel files.

132 xlsx.from.Rdcm

Examples

xlsx.from.Rdcm

Converting .Rdcm files to .xlsx files

Description

A *.Rdcm file contains the list of contents, in dataframe form, of the DICOM files of the same object. The xlsx.from.Rdcm function creates, from a *.Rdcm file, an Excel file, in which each page contains the dataframe representation of a DICOM file of the same object.

Usage

```
xlsx.from.Rdcm(
  Rdcm.filenames,
  dest.dirname = dirname(Rdcm.filenames),
  txt.sep = "\\",
  txt.length = 100,
  tag.dictionary = dicom.tag.dictionary()
)
```

${\tt Rdcm.filenames}$	String vector, representing the *.Rdcm filenames to be converted.
dest.dirname	String vector of the same length as Rdcm. filenames, indicating the directory where the $*.xlsx$ files will be created.
txt.sep	String. Used if as.txt = TRUE. Separator of the tag value elements.
txt.length	Positive integer. Used if as.txt = TRUE. Maximum number of letters in the representation of the TAG value.
tag.dictionary	Dataframe, by default equal to dicom.tag.dictionary, whose structure it must keep. This dataframe is used to parse DICOM files.

xlsx.from.Rdcm 133

Value

Returns a boolean vector, establishing the existence of the created Excel files which have the same basenames as the *.Rdcm files.

```
# First, create a Rdcm file from toy.dicom.raw () to a temporary file for testing.
pat.dir <- file.path (tempdir(), "PM_Rdcm")</pre>
dir.create (pat.dir, recursive = TRUE)
dcm.filename <- tempfile (pattern = "PM_rtplan", tmpdir = pat.dir, fileext = ".dcm")</pre>
zz <- file (dcm.filename, "wb")</pre>
writeBin (toy.dicom.raw (), zz, size = 1)
close (zz)
dicom.to.Rdcm.converter (dcm.filename, pat.dir, update = TRUE)
file.remove (dcm.filename)
list.files (pat.dir)
# Creating an Excel file
Rdcm.filenames <- list.files (pat.dir, pattern = "[.]Rdcm$",</pre>
                              recursive = TRUE, full.names = TRUE)
xlsx.from.Rdcm (Rdcm.filenames)
list.files (pat.dir)
# Cleaning temporary directory
unlink (pat.dir, recursive = TRUE)
```

Index

add.margin, 4, 5, 7, 9, 14, 93	103, 105, 118, 120, 122, 123,
	125–127, 129, 131
bin.closing, 5, 7, 9, 14	
bin.clustering, 6	fan.beam, 55, 56, 58, 59
bin.dilation, 5, 7, 9, 14	fan.planar, <i>55</i> , <i>56</i> , <i>58</i> , <i>59</i>
bin.erosion, 5, 7, 8, 14	fan.sphere, 55, 56, 57, 59
bin.from.roi, 9	fan.to.voxel, <i>55</i> , <i>56</i> , <i>58</i> , <i>58</i>
bin.from.vol, 10, 11	
bin.intersection, 12	get.extreme.pt,59
bin.inversion, 13	get.ijk.from.index,60,67
bin.opening, 5, 7, 9, 14	get.ijk.from.xyz,61
bin.subtraction, 15	get.line,62
bin.sum, 16	get.obj.connection, 43 , 63
castlow str 17 17	get.plane,64
castlow.str, 17, 17	get.rigid.M,66
castup.str, <i>17</i> , 17	get.roi.connection,66,119
dicom.browser, 18, 22, 23	get.value.from.ijk,60,61,67
dicom.parser, 19, 23, 27	get.value.from.xyz,68
dicom.raw.data.anonymizer, 20	get.volume.from.bin, 69, 70
dicom.raw.data.loader, 19, 20, 21	get.volume.from.roi, 69,70
dicom.set.tag.value,22	get.xyz.from.index, 68, 71
dicom.tag.dictionary, 18, 19, 21, 22, 23,	grid.equal,71
25, 26, 80, 81, 83, 131, 132	
dicom.tag.parser, 19, 20, 22, 24	hist, <i>76</i> , <i>77</i> , <i>79</i>
dicom.to.Rdcm.converter, 25, 82, 84, 85,	histo.2D, 28, 50, 54, 72
116	histo.DVH, <i>50</i> , <i>54</i> , 74
dicom.viewer, 26	histo.from.bin, <i>54</i> , <i>75</i> , <i>75</i> , <i>78</i> , <i>80</i>
display. 2D. histo, 27, 73	histo.from.roi, 50, 54, 75, 76, 77, 80
display.3D.contour, 29	histo.vol, <i>54</i> , <i>75</i> , <i>76</i> , <i>78</i> , <i>79</i>
display.3D.mesh, 30	
display. 3D. sections, 31	load.obj.data, 80, 81, 82, 84, 85
display.3D.stack, 32	load.obj.from.dicom, 80, 81, 82, 84, 85
display.dV_dx, 36, 38, 76, 78, 80	load.obj.from.Rdcm, 48, 54, 80, 81, 82, 84,
display.DVH, 34, 36, 75	85
display.DVH.pc, 34, 35, 75	load.patient.from.dicom, 42, 54, 63, 66,
display.histo, 37, 37, 76, 78, 80	68, 77, 83, 85, 89, 94, 98, 99, 101,
display.kplane, 39, 47, 60	115, 120, 128, 129
display.legend, 41	load.patient.from.Rdcm, 10, 30, 32, 33, 42,
display.obj.links, 42, 64	46, 54, 60, 63, 66, 68, 77, 84, 84, 89,
display.palette, 43	94, 98, 99, 101, 106, 109, 115, 118,
display.plane, 39, 40, 45	120, 128, 129
	load.Rdcm.raw.data,85
espadon.class, 4-8, 10, 12-16, 28-30, 32,	load.T.MAT, 10, 30, 32, 33, 46, 48, 54, 60, 66,
<i>35</i> , <i>37</i> , <i>39</i> , 48, <i>55</i> , <i>56</i> , <i>58</i> , <i>81–84</i> ,	68, 77, 84, 85, 86, 89, 94, 98, 99,

INDEX 135

```
101, 106, 109, 115, 118, 120, 128,
         129
mesh.from.bin, 30, 31, 50, 54, 87
mesh.in.new.ref, 89
mesh.repair, 90
mesh.spheric.proj, 90
nesting.bin, 92
nesting.cube, 4, 5, 7, 9, 14, 93
nesting.roi, 94
orientation.create, 95
pal.RVV, 96
plot, 34, 36-38
Rdcm.inventory, 97
ref.add, 89, 98, 99-101, 115, 120, 128, 129
ref.cutplane.add, 98, 99, 101, 115
ref.remove, 98, 99, 100, 101
ref.srctodest.add, 101
rt.chi.index, 102, 105
rt.gamma.index, 103, 104
rt.indices.from.bin, 106, 114
rt.indices.from.roi, 107, 108
save.T.MAT, 115
save.to.Rdcm, 116
select.names, 70, 117
shade3d, 30
struct.clustering, 66, 67, 118
struct.from.bin, 119, 122
struct.in.new.ref, 120, 128
struct.merge, 121
toy.dicom.raw, 122
toy.load.patient, 42, 54, 63, 123
vcgClean, 87
vcgSmooth, 87, 88
vector.product, 124
vol.copy, 124
vol.create, 125
vol.from.bin, 126
vol.gradient, 127
vol.in.new.ref, 121, 128
vol.median, 128
vol.regrid, 129
vol.sum, 130
xlsx.from.dcm, 27, 131
xlsx.from.Rdcm, 132
```